Product Information AlproJet-W

Cleaning of dental aspiration systems

Highly efficient aldehyde-free and nonfoaming liquid concentrate for the weekly contra-cleaning and maintenance of dental aspiration and separation systems, with and without amalgam separator.

- supplementary product with special protein dissolvent maintaining the full aspiration capacity of the system, even in case of high protein load, e.g. before and after surgical treatments
- should be used at least two days per week (e.g. Monday and Tuesday) in order to split, dissolve and flush away any residues
- dissolves lime (acid pH value), neutralizes hydrogen peroxide H₂O₂ and sodium bicarbonate (from airflow)
- free from phenols, aldehydes, phosphates and chlorine!
- according to scientific examinations AlproJet-W releases even less mercury from the separated amalgam than tap water

Composition

Amidosulfonic acid, tensides, auxiliary substances

Efficacy

Cleaning

Concentrations

5 % with warm water ca. 40° - 50° C (in case of heavy contamination up to 10 %)

Clinical data

 Prof. Dr. P. Heeg, Tübingen, screening test, 1996-08

Please feel free to request these clinical data from us.

Field of application

Liquid concentrate for the weekly cleaning and maintenance of aspiration and separation systems, with and without amalgam separator.

MANUFACTURER: ALPRO MEDICAL GMBH Mooswiesenstraße 9 78112 St. Georgen Germany



as at: 15 July 2015





Shelf life

3 years

Application

For the reliable cleaning and maintenance it is recommended to use AlproJet-W 2 days a week instead of AlproJet-D or AlproJet-DD (ARTG No. 164974), in each case 1 – 2 times a day (noon, evening). In case of high protein load a daily application of AlproJet-W is also temporarily possible.

2 litres of a 5 % ready for use solution are required for one application. For this purpose dose 2 x 50 mL AlproJet-W concentrate by using the dosage chamber of the dosage bottle and fill up in the AlproJet MixCup with $40 - 50^{\circ}$ C warm water to 2 litres of a ready for use solution. Afterwards aspirate 1 litre of the ready for use solution through the aspiration tubes and pour the remaining 1 litre into the spittoon bowl.

Thereafter let take effect for at least 30 minutes and rinse the aspiration tubes and the spittoon bowl thoroughly with water just before starting work.

Application only by qualified personnel wearing protective clothing.

2 +49 7725 9392-0
⇒ +49 7725 9392-91

Product Information AlproJet-W



Cleaning of dental aspiration systems

as at: 15 July 2015

Environment

AlproJet-W does neither contain chlorine nor oxygen splitters and is free from aldehydes and phenols. Biodegradable according to OECD regulations.

Approvals

Approved and recommended amongst others by Sirona Dental.

Material compatibility

AlproJet-W is highly compatible with materials such as plastics, rubber and non-ferrous metal, alloyed steel etc.

Due to the excellent material compatibility AlproJet-W may be used in the dental unit overnight and/or over the weekend.

SAFETY DIRECTIONS

CORROSIVE - Contains sulfamic acid 4-5 % w/w

- Avoid contact with eyes
- Avoid contact skin

Delivery forms

REF 3106	1 L Dosage bottle
REF 3105	5 L Can

Accessories

er tap
pump
nch
pump for 1 ltr. bottle
MixCup

Distributor:

Biodegree Pty Ltd Sydney, Australia