Science.
Applied to Life.™

STEP-BY-STEP



3M™ RelyX™ Universal





System Components

3M™ RelyX™ Universal

Resin Cement

• Dual-cure two paste resin cement

Adhesive and self-adhesive indications




3M™ Scotchbond™ Universal Plus

Adhesive

- Self-etch, selective-etch and total-etch adhesive
- MDP primer for zirconia and metal
- Silane primer for glass ceramics, composites, fiber reinforced posts

Shade Selection



3M™ RelyX™ Universal Micro Mixing Tips with Elongation Tips









Self-sealing automix syringe







Full rotation until stop.



Apply moderate centric force without bending the plunger.



Discard a small amount of cement.







Store syringe without mixing tip.



Restoration Pre-Treatment Overview

Fiber-reinforced Composite Post

Priming with adhesive not required for $3M^{\mathsf{TM}} \operatorname{RelyX^{\mathsf{TM}}}$ Fiber Posts





Zirconia, Alumina, Metal

Priming with adhesive optional, recommended for veneer, table top, adhesive bridge





See section

1

Composite/Hybrids







Glass-Ceramic

Alternative primer: Apply a Silane, e.g. 3M™ RelyX™ Ceramic Primer









Follow the instructions for use of the restoration material for specific instructions. **Option**

Self-adhesive. Recommended for:

- Post
- Crown
- Bridge



Adhesive on tooth not required.

See section

2a/2b



(Selective-etch) adhesive. Recommended for:

- Inlay
- Onlay





See section

Total-etch adhesive. To be used for:

- Table Top
- Veneer
- Adhesive **Bridge**





See section


Restoration Pre-Treatment Recommendation



Fiberreinforced Composite Post







Zirconia, Alumina, Metal









Y-STEP

3M™ RelyX™ Universal Resin Cement

Composite/ Hybrids











Glass-Ceramic












Self-Adhesive Post Cementation

STEP-B













To unlock syringe fully rotate mixer until stop.






Do not use high pressure. Check mixer position (image 5) if extrusion force is high.



No lentulo!















Self-Adhesive Cementation

STEP-B

1

Pre-treat restoration according to section 1









To unlock syringe fully rotate mixer until stop.



6

Do not use high pressure. Check mixer position (image 5) if extrusion force is high.









Tack cure approx. 2-3 sec per surface (>1200 mW/cm²).

10





Light cure 10 sec per surface (>1200 mW/cm²) or wait 6 minutes.



13



(Selective-Etch) Adhesive Cementation

STEP-B

1

Pre-treat restoration according to section 1




Selective-etch of uncut enamel. Optional for cut enamel.











To unlock syringe fully rotate mixer until stop.



Do not use high pressure. Check mixer position (image 7) if extrusion force is high.









Tack cure approx. 2-3 sec per surface (>1200 mW/cm²).



Light cure 10 sec per surface (>1200 mW/cm²) or wait 6 minutes.







Total-Etch Adhesive Cementation

STEP-B

1

Pre-treat restoration according to section 1



Total-etch.













To unlock syringe fully rotate mixer until stop.



Do not use high pressure. Check mixer position (image 7) if extrusion force is high.

9






11



Tack cure approx. 2-3 sec per surface (>1200 mW/cm²).

12



Light cure 10 sec per surface (>1200 mW/cm²) or wait 6 minutes.



14





Please refer to product instructions for use for complete information.



3M Australia Pty Limited Building A, 1 Rivett Rd North Ryde, NSW 2113 1300 363 454 www.3M.com.au/RelyXUniversal 3M New Zealand Limited 94 Apollo Drive Rosedale, Auckland 0632 0800 808 182 www.3M.co.nz/RelyXUniversal

Always follow the Instructions for Use (IFU) and refer to IFU for full indications, precautions and warnings. 3M, "3M Science. Applied to Life.", RelyX and Scotchbond are trademarks of 3M. © 3M 2020. All rights reserved.