


# NINJA II

## Technical specifications

Intended use	Odontotomy			
Indications for use	Root cutting			
Description	Saw with double grinding, marked every 3 mm, for hemisections and radicular amputations			
Tip type	Saw			
Allergenic factor	Tip made of surgical stainless steel as per ASTM F899 - 12b  Surgical stainless steel contains: <ul style="list-style-type: none"> <li>• nickel</li> <li>• chromium</li> </ul>			
Commercial reference	NINJA II : F87545			
Medical class	Class of medical device: IIa according to 93/42/EEC directive			
Adjustments	Ultrasonic generators	Maximum mode	Maximum fine setting	Irrigation
	Piezotome Cube	D1	---	60-80 ml/min
	Implant Center 2 Piezotome Solo LED	D1	1-3	60-80 ml/min

## Service life

As it is not possible to establish a maximum number of uses (that may be determined by many parameters such as duration of use, hardness of tissue, the force applied, wear), we recommend that medical devices being routinely used are replaced at least once a year. Replace the medical device if the tip does not vibrate at the expected frequency, the treatment is not progressing as normal and is taking longer or is at a standstill.

## Associated documentation


Document title	References
Consulting electronic user instructions	J00007
Cleaning, disinfection and sterilisation instructions for tips	J02001
Extraction Kit clinical tips	J57351
Dental surgery ultrasonic generator power settings table	J58010
General instructions relating to the complete range of dental surgery tips	J87301

# NINJA II


## Sterilisation




Sterilise prior to each use




The tips must be cleaned, disinfected and sterilised prior to each use. Prior to each use, please refer to the cleaning, disinfection and sterilisation instructions for tips J02001.



## Installation



## Disassembly



## Manufacturer identification



SATELEC  
A Company of ACTEON Group  
17, avenue Gustave Eiffel  
ZI du Phare  
33700 MERIGNAC  
France  
Tel. +33 (0) 556.34.06.07  
Fax. +33 (0) 556.34.92.92  
E.mail : satelec@acteongroup.com  
www.acteongroup.com



## Date of first CE marking

2009



The practitioner has to continually monitor the clinical site to immediately identify any heating risk.

Tip | NINJA II | J01011 | (09) | V5 | 05/2021 | NO50EN060E

SATELEC | A Company of ACTEON Group  
17 av. Gustave Eiffel | ZI du Phare  
33700 MERIGNAC | FRANCE  
Tel. +33 (0) 556 34 06 07 | Fax. +33 (0) 556 34 92 92  
satelec@acteongroup.com | www.acteongroup.com

Download



Instructions For Use

