EcuSphere

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English

Product description

The EcuSphere System consists of four light-curing composites specifically developed for their respective applications:

The radiopaque, high viscosity EcoSphere-Carat is preferred for permanent restorations, especially in the posterior area.

The radiopaque EcuSphere-Shape is best suited for modeling, as dentine material for high quality aesthetic anterior tooth care and as a universal composite.

EcuSphere-Shine is preferred for use as an artificial enamel in the anterior tooth area.

Radiopaque EcuSphere-Flow is suitable for all indications that require an easy-flowing composite and it perfectly matches the other EcuSphere colors.

The color coding of the individual materials on the syringe shaft or Safetip plunger ensures easy distinction (EcuSphere-Carat = black, EcuSphere-Shape = grey, EcuSphere-Shine = ivory, EcuSphere-Flow in the Luer-Lock syringe).

Indications

EcuSphere-Carat	Restorations of all cavity classes I-V, production of direct inlays and onlays
EcuSphere-Shape	Restorations of all cavity classes I-V
EcuSphere-Shine	Restorations in the anterior tooth area of cavity classes III, IV and V, production of direct veneers
EcuSphere-Flow	Small restorations of cavity classes I, III, IV, V (cervical caries, root erosions, wedge-shaped defects), milk tooth cavities, minimally invasive fillings, liners, expanded fissure sealing, blocking out and filling of undercuts in cavities, repairs of composite restorations

Contra-indications

- ▶ Application of EcuSphere is contra-indicated if a dry working area or the recommended application technique is not possible. Do not use directly on exposed pulp or in cases of allergies to any of the ingredients.

Basic safety instructions

- ▶ For dental use only!
- ▶ Keep out of reach of children!
- ▶ Avoid contact with the skin! In the event of inadvertent skin contact, wash the affected area thoroughly with soap and water. Wear suitable protective gloves.
- ▶ Avoid contact with the eyes! In the event of inadvertent contact with the eyes, rinse thoroughly with plenty of water immediately and consult a physician if necessary.
- ▶ Avoid contact with oral soft tissue. In the event of inadvertent contact, remove the material immediately and rinse the affected area with plenty of water. If irritation of the soft tissue persists, consult a physician.

Side effects

There have been no systemic side effects to date. Isolated cases of hypersensitivity or contact allergy to components of the material cannot be excluded.

Interactions with other materials

- ▶ Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil.
- ▶ The use of cationic mouthwashes, plaque indicators and chlorhexidine may lead to discoloration.

Instructions for use

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm². Position the curing light as close to the material as possible.
- ▶ Do not remove the oxygen inhibited layer because this is required for the bond with the next layer.
- ▶ Observe the manufacturer's instructions for any other materials that are used with EcuSphere.

Recommended use

Preparation

1. Before preparation, clean the tooth and remove all cleaning residue with water spray. Determine VITA tooth color.

Note: DMG recommends the use of a rubber dam.

2. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Protect the areas close to the pulp with a suitable material.
3. Use matrices to protect the neighbouring teeth.
4. Only when using the Total-Etch technique: Etch the preparation in accordance with the manufacturer's instructions. DMG recommends TECO as an adhesive. There is no etching step when using the Self-Etch technique. In this case, DMG recommends the self-etching adhesive Contax.
5. Apply bonding agent in accordance with the manufacturer's instructions.

Application

Note: Apply all four EcuSphere materials in layers of a maximum of 2 mm thick each and light cure each layer for 20 s.

1. If necessary, apply a liner of thin-flowing composite or compomer (e.g. EcuSphere-Flow or Ionosit-Baseliner).
2. Fill the cavity with EcuSphere-Carat or EcuSphere-Shape in the desired color, in layers of a maximum thickness of 2 mm. Light cure each layer for 20 s or mold dentine core out of opaque material, e.g. EcuSphere-Shape or EcuSphere-Carat in the desired shade.
3. For highly aesthetic restorations in the anterior tooth area, simulate the enamel with EcuSphere-Shine.

Processing

- ▶ Remove excesses, shape and finish, e.g. with diamond finishing burs and flexible disks with reducing grit sizes.
- ▶ Carry out high shine polishing of EcuSphere-Carat, EcuSphere-Shape and EcuSphere-Shine with polishing pastes using polishing cups/disks.

Disinfection of the Safetip-Dispenser

The Safetip-Dispenser is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfecting solutions listed by the RKI (Robert Koch Institute). The Safetip-Dispenser is also autoclavable.

Storage and shelf life

- ▶ Store in a dry place at room temperature (15 to 25 °C/59 to 77 °F)!
- ▶ Reseal syringes after use!
- ▶ Do not use after the expiration date!

Composition

EcuSphere-Carat

Dental glass in a Bis-GMA based matrix of dental resins, pigments, additives, catalysts. Filler content: 77 % by weight = 57 % by vol. (0.02 – 0.7 µm)

EcuSphere-Shape

Dental glass in a Bis-GMA based matrix of dental resins, pigments, additives, catalysts. Filler content: 77 % by weight = 57 % by vol. (0.02 – 1.5 µm)

EcuSphere-Shine

Micro-filler and micro-filled prepolymers in a Bis-GMA based matrix of dental resins, pigments, additives, catalysts. Filler content: 76 % by weight = 55 % by vol (0.04 – 20 µm).

EcuSphere-Flow

Dental glass in an optimised matrix of Bis-GMA; 63 % by weight = 41 % by volume inorganic filler (0.02 – 3 µm).

Classification

ISO 4049:2009

Curing depth (20 s) for EcuSphere-Carat: 1.75 mm.

Curing depth (20 s) for EcuSphere-Shape: 1.75 mm, opaque colours = 1.5.

Curing depth (20 s) for EcuSphere-Shine: 1.75 mm.

Curing depth (20 s) for EcuSphere-Flow: 1.75 mm, opaque colours = 1.25 mm.

Packaging

Starter Kit	
1 syringe containing 3 g EcuSphere-Carat (A3.5) 1 syringe containing 3 g EcuSphere-Shape (D-A3*) 1 syringe containing 2.5 g EcuSphere-Shine (Transparent) 1 syringe containing 1.5 g EcuSphere-Flow (A3.5)	REF 220324

EcuSphere-Carat

Packaging		
1 syringe containing 3 g paste	A1 = REF 220287 A2 = REF 220267	A3 = REF 220285 A3.5 = REF 220266
20 Safetips each containing 0.3 g paste	A1 = REF 220323 A2 = REF 220271	A3 = REF 220322 A3.5 = REF 220270

EcuSphere-Shape

Packaging		
1 syringe containing 3 g paste	A1 = REF 220334 A2 = REF 220308 A3 = REF 220299 A3.5 = REF 220297 A4 = REF 220347 B1 = REF 220335	B2 = REF 220336 B3 = REF 220337 C1 = REF 220338 C3 = REF 220339 D2 = REF 220340 D-A3* = REF 220291
20 Safetips each containing 0.3 g paste	A1 = REF 220327 A2 = REF 220309 A3 = REF 220307 A3.5 = REF 220298 A4 = REF 220346 B1 = REF 220328	B2 = REF 220329 B3 = REF 220330 C1 = REF 220331 C3 = REF 220332 D2 = REF 220333 D-A3* = REF 220292

EcuSphere-Shine

Packaging		
1 syringe containing 2.5 g paste	Transparent = REF 220262	A2 = REF 220261 A3 = REF 220260
20 Safetips each containing 0.3 g paste	Transparent = REF 220258	A2 = REF 220257 A3 = REF 220256

EcuSphere-Flow

Packaging		
2 syringes each containing 1.5 g, 10 Luer-Lock-Tips	A1 = REF 220315 A2 = REF 220314	A3 = REF 220313 A3.5 = REF 220312

Accessories

60 Luer-Lock-Tips, black	REF 211477
1 mixing block, 62 x 62 mm	REF 212988
1 Safetip-Dispenser	REF 213122

* D-A3 = Dentine opaque A3

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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