



# COLGATE® PULSE

## SERIES 1

Tracks how long and often you brush to improve your brushing routine

2 modes, 2 minute timer and 2 week battery life†



CG-61030574

CG-61030573

## SERIES 2

Tracks how long, often & well you brush, coaching you to a healthier smile

3 modes, pressure sensor, 2 minute timer, 2 week battery life† and travel case



CG-61030577

CG-61030576



## REPLACEMENT BRUSH HEADS

Not available at Henry Schein. Only available across select national retailers.



**SENSITIVE**

Tough on plaque, gentle on gums



**WHITENING**

Whiter teeth in 1 week\*



**DEEP CLEAN**

Removes 5x more plaque†



# PULSE

SERIES 2

## SMART COACHING FOR A SUPERIOR CLEAN\*



Available from



\*at gum margin after brushing for 6 weeks with vs without Colgate Connect app



Available from



Visit [henryschein.com.au](https://henryschein.com.au) for more details


†when used twice daily  
\*vs. an ordinary manual toothbrush  
†at the gum margin after 4 weeks vs a manual flat-trim toothbrush

MOST PEOPLE MISS OVER **40%** OF THEIR TEETH WHEN THEY BRUSH\*


GET CONNECTED AND BRUSH BETTER



**GUIDED BRUSHING**  
16 zone detection with speed and angle feedback (Series 2 only)



**SMILE REWARDS**  
Get into a healthy brushing routine and earn points#



**GET REMINDED**  
Engaging daily in-app content and reminder notifications

## FEATURES



**3 BRUSHING MODES**  
A simple switch lets you select the right mode for you: Gentle Clean, Daily Clean, Deep Clean (Series 2 only)



**CELEBRATE BRUSHING STREAKS**  
The Smile Light changes colours to indicate brushing achievements and to encourage more consistency



**PRESSURE SENSOR**  
Uses smart technology to detect over-brushing, red LED lights alert you when brushing too hard (Series 2 only)



**2 MINUTE TIMER**  
White LED lights help guide your pace for a complete 2 minute brushing routine



**RECHARGEABLE BATTERY**  
2 week battery life when used twice daily

REMOVE **50% MORE PLAQUE\*\*** (SERIES 2 ONLY)

\*Based on tooth surface coverage data collected by the Colgate Connect app when brushing with a smart electric toothbrush.  
\*\*at gum margin after brushing for 6 weeks with vs without Colgate Connect app.

#See colgate.com/apprewards for more information on how you could get rewarded