



roeko

Surgitip

Aspirator Tips

Meet your aspirations in surgery
and endodontics

 COLTENE




roeko

Surgitip

Surgitip aspirator tips are proven aids in drying the oral field of treatment in surgery and endodontics. They are intended for low volume aspiration of liquid and solid substances and enable aspiration in restricted dental operating environment. Their specific design make them suitable for working under microscope and loupe. Hygienic handling properties and ease of use make them an enjoyable device in the dental office.

- Green tip provides excellent contrast to the treatment site
- Open field of view due to the narrow shape and conical tip
- Hygienic and time-saving, being sterile, individually packaged and ready for immediate use
- Various tip geometries and diameters for a wide range of treatments



Surgitip-endo

- 0.35 mm tip opening
- sterile

Surgitip-endo is an aspirating tip intended for aspiration during endodontic treatments and is specific for drying of the root canal. Surgitip-endo removes rinsing solutions and moisture in the root canal. Less paper points are then needed for drying the root canal. The specific ball design of the suction tip allows for a multi-directional flexibility. It can easily be introduced into hard-to-reach root canals without having to bend the Canal Tip. Aspiration is highly effective at all angles and visibility is not constricted.

Surgitip-micro


- 1.2 mm tip opening
- sterile

Surgitip-micro is intended for low volume aspiration of liquid and solid substances during dental surgical and micro-surgical procedures in the oral cavity. The cannula offers good grip and is easy to manoeuvre which enables targeted aspiration in restricted dental operating environments.



Surgitip

- 2.8 or 4.9 mm tip opening
- sterile



Surgitips are intended for low volume aspiration during dental surgical procedures. With either 2.8 or 4.9 mm it can be used whenever greater aspiration capabilities are required. The 4.9 mm tip consists of the bent aspirator cannula without the green tip.

Aspirator tips


- 4.9 mm tip opening
- non sterile

Aspirator tips are non sterile and intended for low volume aspiration of liquid and solid substances during non-surgical dental procedures. They come in bulk packs with 100 pieces.




Double Adapter

- end opening 11 and 16 mm
- reprocessable



Double Adapters connect aspiration cannulas to standard suction systems with ø 11 mm or ø 16 mm openings.




ORDER INFORMATION

	Tip opening \varnothing mm	End opening \varnothing mm	Contents	REF
Surgitip-endo	0.35 mm	6.5 mm	20 pcs *	60018280
Surgitip-micro	1.2 mm	6.5 mm	20 pcs *	462012
			100 pcs *	462013
Surgitip	2.8 mm	6.5 mm	20 pcs *	462025
			100 pcs *	462026
			20 pcs *	462048
	4.9 mm	6.5 mm	100 pcs *	462049
Aspirator tips	4.9 mm	6.5 mm	100 pcs *	461048
Adapter		11 mm	1	463025
		16 mm	1	463011
Double Adapter		11 + 16 mm	1 Set	463000

* incl. 1 Double Adapter

© COLTENE – www.coltene.com



Coltene/Whaledent Ltd.
The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF/GB
T +44 1444 235486
info.uk@coltene.com

Coltene/Whaledent AG
Feldwiesenstrasse 20
CH-9450 Altstätten
T +41 71 757 5300
info.ch@coltene.com

Coltene/Whaledent GmbH + Co. KG
Raiffeisenstraße 30
DE-89129 Langenau
T +49 7345 805 0
info.de@coltene.com