

## Powder-Free Nitrile Examination Gloves

### Strong and puncture resistant

- Latex free, Powder free
- Excellent tactile sensitivity, finger textured
- 240mm length with beaded cuff
- Ambidextrous, single use, non-sterile
- Conforms to EN 455-1,2,3,4. ARTG # 320857



### Technical Data:

Length	≥ 240mm	
Thickness (@palm)	≥ 0.06 mm	
Colour	Blue	
Force at Break (N)		
Unaged	≥ 6.0	
After ageing	≥ 6.0	
Watertight test	AQL 1.5	Inspection level G1
Visual Defects Major	AQL 2.5	Inspection Level G1
Powder residue	Max 2.0 mg/glove	

### Complies with

Medical Device Regulation 2017/745. EN 455-1,-2,-3,-4. EN 420:2003 +A1

Personal Protective Equipment Regulation (EU) 2016/425 -Category 1, for minimal risk.

### Packaging:

200 gloves per box. Carton of 10 boxes.

9796093	XSmall	200 per box
9796094	Small	200 per box
9796095	Med	200 per box
9796096	Large	200 per box

### Storage:

Store below 30° C in a cool, dry area away from direct sunlight.