

# MAXI-GARD PHOSPHOR BARRIER ENVELOPES



Barrier envelopes are a disposable, one use item which guarantees effective hygiene when handling a phosphor plate and placing them in the patient's mouth.

The Maxi-Gard Barrier Envelopes are an economical solution to ensure that the phosphor plate is not in contact with saliva and/or the patient.

## Maxi-Gard VL (Vertical Loading) Phosphor Plate Barrier Envelopes

- Accepts and unloads plates vertically.
- Perfect for PSP machines that accept plates along the short side.
- Plates load quickly and easily.
- Butterfly seam makes plate removal fast and easy.
- Semitranslucent on one side so you can verify proper plate position.

Compatible with units that accept **vertical loading** such as:

Durr Vistascan range	Carestream 7600
Air Techniques Scan X duo	Sirona Xios
Acteon PSPIX	My Ray Hy Scan

HS-5701558	Maxi-Gard VL Barrier Envelope	Size 0	Box 100
HS-5701559	Maxi-Gard VL Barrier Envelope	Size 1	Box 100
HS-5701560	Maxi-Gard VL Barrier Envelope	Size 2	Box 300



## Maxi-Gard HL (Horizontal Loading) Phosphor Plate Barrier Envelopes

- Accepts and unloads plates horizontally.
- Plates load quickly and easily.
- Black-on-black material for maximum protection.
- Excellent quality at a great value.

Compatible with units that accept **horizontal loading** such as:

Air Techniques Scan X duo	Sirona Xios
Carestream 7200	My Ray Hy Scan
Apixia PSP	

HS-5701561	Maxi-Gard HL Barrier Envelope	Size 0	Box 100
HS-5701563	Maxi-Gard HL Barrier Envelope	Size 1	Box 100
HS-5701564	Maxi-Gard HL Barrier Envelope	Size 2	Box 300



# MAXI-GARD PHOSPHOR BARRIER ENVELOPES



## Maxi-Gard MC (Magnet/Cardboard) Phosphor Plate Barrier Envelopes

- For PSP plated with magnets.
- Accepts and unloads plates vertically.
- Comes with cardboard inserts to protect plates.
- Butterfly seam for effortless plate removal.
- Semitranslucent on one side so you can verify proper plate position.

Compatible with units that require **cardboard** protection and plates that have **magnets** such as:

Acteon PSPIX	Gendex
Kavo Scan Exam	Soredex Optime/Digora
Instrumentarium Origo	Sirona Xios



HS-5701565	Maxi-Gard MC Barrier Envelope	Size 0	Box 200
HS-5701566	Maxi-Gard MC Barrier Envelope	Size 1	Box 200
HS-5701567	Maxi-Gard MC Barrier Envelope	Size 2	Box 200