

ECOLAB[®] MICROFIBRE TOWEL

Lint Free Knitted Polyester / Polyamide

Washable | Autoclavable | Super-absorbent

40cm x 40cm | Blue or Yellow

Pack of 5

EC-EMTB5BLUEE | EC-EMTB5YELLOWE



EcoLab Microfibre towels are perfect for sterilising rooms in the drying/packaging area and in the surgery.

The polyester towels are designed to absorb excess moisture and streamline the sorting process.

With two colours it's easy to identify the different zones within the packaging area:

Blue – can be identified as the drying towel

Yellow – can be identified as the packaging towel

By ensuring the segregation between colours, this assists staff with the workflow of the sterilisation room and helps them remember to remove wet gloves before moving to the packing towel so clean and dry instruments can be packaged appropriately.

COMMON USES INCLUDE

- Bench Top Towels for Instruments
- Washer for Dental Chairs to remove Detergent Residue
- Wipe for Cleaning Glass or Lensed Sensitive Instruments
- General Cleaning Purposes



WASHING INSTRUCTIONS

Preferred - Washing machine cycle (up to 95 degrees)

Recommended - Cleaning with a detergent solution, rinsed in hot water and hung to dry overnight

*Can be autoclaved but this will not remove dirt – will require pre cleaning. If the intended purpose is to reuse in packaging area, then the above washing methods are preferred over autoclaving.