

DENTIQ

Your Chairside 3D Printer



Engineered to endure the most demanding dental workflows. Precision, versatility and, reliability for the Dentistry of the Future. Enabled with the latest LCD technology and UV backlight to deliver high quality results consistently.



Ideal for **same-day dentistry**



Auto-calibration and **100+ Prevalidated Materials**



Intuitive and **easy to use**



47 μ m Accuracy and consistency



Premium Ownership Experience after purchase

TECHNICAL SPECIFICATIONS

3D Print Technology	LCD
Machine Size	25 x 23 x 38 cm
Weight	10.5 Kg
Print Size	120 x 68 x 140 mm
XY Resolution	47 μ m
Slice Thickness	30 - 150 μ m
Resin Wavelength	385 - 405 nm
Power Supply	110 - 240 V AC, 50 - 60 Hz

PACKAGE INCLUDES

- 2 x Extra Print Screens
- 1 x Extra Vat
- 1 x Ethernet Cable
- 1 x Wi-Fi Dongle
- 2-Year Warranty



FEATURES

- Precalibrated 100+ resins
- Smaller footprint ideal for clinics
- 47µm pixels for smoother surfaces
- Auto-calibration of Z Axis
- Replace most parts in under 15 min.
- Dental focused technical support
- Intuitive UI for quick implementation
- 2-Year Warranty



APPLICATIONS



Dental Models
2 per print
100 min



Surgical Guides
9 per print
90 min



Temporary Crowns
40 per print
50 min



Denture Bases
2 per print
140 min



Denture Teeth
2 per print
60 min



Removable Partial Dentures
3 per print
100 min



Gingiva
32 per print
75 min



Impression Trays
1 per print
160 min



IBT
3 per print
210 min



Castable Crowns
40 per print
50 min



Mouth Guards
2 per print
160 min



Splints
2 per print
160 min

ACKURETTA ECOSYSTEM

Envision a future where dental clinics have a solution to 3D print end-product restorations chairside with lower cost, and higher efficiency. From Hardware, Software, Finishing Solutions, Materials, to Training Resource, **Ackuretta keeps everything In-house and completed in the same day.**

Software



3D Printing Materials



Finishing Kit



Industrial-grade 3D printers

