

## DEHP Gloves Latex Powder-Free

**Non-sterile powder-free latex gloves. Textured fingers.**

- Finger textured for enhanced grip
- 240mm length with beaded cuff
- Chlorinated to make donning easier
- Ambidextrous, single use, non sterile
- Complies with AS/NZS 4011.1:2014, ARTG # 320857



### Technical Data:

Tested as per requirements of AS/NZS 4011:2014

Length:  $\geq 240\text{mm}$

Thickness (@palm)  $\geq 0.08\text{mm}$

Weight: 5gm (Size M)

Force at Break

Before ageing  $\geq 7.0\text{N}$

After ageing  $\geq 6.0\text{N}$

Total Proteins  $\leq 50 \mu\text{g/g}$

Watertight test. AQL 1.5 Inspection level G1

This product contains natural latex and dithiocarbamate and can cause allergic reaction. Do not use if sensitive to these substances.

### Packaging:

Cartons of 10 boxes. 200 gloves per box.

HSD-9770361	XSmall	200 per box
HSD-9770362	Small	200 per box
HSD-9770363	Med	200 per box
HSD-9770364	Large	200 per box

### Storage:

Store below 30° C in a cool, dry area away from direct sunlight.