

WATERLINE DISINFECTION



**HIGH QUALITY WATERLINE PROTECTION
AT GREAT VALUE**

NEW Henry Schein® Dental Unit Water Purification Cartridges and Waterline Maintenance Tablets

KEEP PRACTICES SAFE AND COMPLIANT WITH CDC AND ADA STANDARDS.

DENTAL UNIT WATERLINE TREATMENT PRODUCTS

Henry Schein® Dental Unit Water Purification Cartridge

2 Cartridges/Box 570-2562

Henry Schein® Waterline Maintenance Tablets (50 Tablets/Box)

2L 570-2560 | 750mL 570-2561

WATERLINE DISINFECTION



HIGH QUALITY WATERLINE PROTECTION

NEW Henry Schein® Waterline Maintenance Tablets help practices exceed the ADA and CDC standards for clean water by 50x.

TRUST OUR EASY-TO-USE TABLETS FOR RELIABLE SAFETY AND COMPLIANCE.

DENTAL UNIT WATERLINE MAINTENANCE TABLETS

ADVANTAGES

- **Compliant** - Tablets exceed the ADA and CDC standards for clean water by 50X, maintaining 10 CFU/mL or less, when used as instructed.
- **Easy to Use** - Simply place a tablet in your empty dental unit waterline bottle every time you refill it.
- **Effective** - Tablets may be left in waterlines overnight and for extended periods of non-use. They have been shown to be effective for up to 28 days.
- **Safe, Odorless, and Tasteless** - Tested for no ill-effect on bonding or equipment.

50 Tablets/Box 2L 570-2560 | 750mL 570-2561

WATERLINE DISINFECTION



HIGH QUALITY WATERLINE PROTECTION

NEW Henry Schein® Dental Unit Water Purification Cartridge two-cartridge system treats waterlines with 40% greater power when compared to single-cartridge systems.

EASILY KEEP YOUR PRACTICE SAFE AND COMPLIANT.

DENTAL UNIT WATER PURIFICATION CARTRIDGE

FEATURED ADVANTAGE

Our dual cartridge system's total active antimicrobial released over the course of 1 year (or 240L) is 40% more than traditional single cartridge systems.

This means practices can achieve improved water test results plus more consistent compliance with CDC and ADA drinking water standards (i.e., ≤ 500 CFU/mL of heterotrophic water bacteria).

ADDITIONAL BENEFITS

- **Complimentary Reusable Dummy Straw** (Never shock through a cartridge. It will ruin it.)
- **Convenient Replacement Indicators** - Kit includes two 6-month automatic replacement indicators that act as a visual reminder for when to replace or reorder.
- **Simple Instructions** - Installs in minutes and can be used with tap or distilled water.
- **Safe for Patient Contact** - Non-toxic, non-allergenic, and non-corrosive.

2 Cartridges/Box 570-2562