

Electrodes



I40S I40CA
INCISION
Ø 0.40 mm (RODS)



I22S I22CA
INCISION
Ø 0.22 mm (NEEDLES)

Abscess incision, exposure of impacted tooth, cutting of peri-coronary sac, hemostasis, exposure of crown lengthening, gingivoplasty, gingivectomy, frenectomy, tissue ablation before impression, stripping, etc.



TR22R TR22L TR22T
EXCISION
Ø 0.22 mm (LOOPS)



FC10N FC25B FC32B
INCISION/
FULGURATION
Ø 1 mm CONICAL)
FULGURATION/COAGULATION
Ø 2.5 mm - Ø 3.2 mm
(SPHERICAL)