




STS

PROFESSIONAL
CREATING PRODUCTS.

HART TEST

HOLLOW AIR REMOVAL TEST FOR SMALL STEAM STERILISERS
TO ISO 11140-6:2022

STS Professional's
Hollow Air Removal Test or
'HART' test (ARTG 460985)
is manufactured to
ISO 11140-6:2022 in the
'helix' style to demonstrate
the air removal & steam
penetration efficiency for
hollow items in small
steam sterilisers.



HART
HOLLOW AIR REMOVAL TEST
**ROUTINE &
DAILY PCD
ISO 11140-6**

The HART test is made up of the
PCD & 250 Type 2 chemical indicator strips
& may be used for **daily testing &/or routine
batch monitoring** in every cycle for hollow
air removal & steam penetration testing to

AS5369:2023.

More details &
instructions for use
found enclosed.



REORDER CODE: 9S-HART
STSPRO.COM.AU