



Traverse™

Rotary Glide Path File

Own the glide path.

Introducing the Traverse rotary glide path file. Over 125 years ago, Kerr invented the K-file and now we unveil the latest, state-of-the-art file technology — Traverse. Simplifying the procedure, while allowing more control, Traverse inspires confidence and safety, providing a tapered glidepath, efficient, minimally invasive procedures and optimum results.

Efficiency + Control = Piece of Mind.



Resistance To Stress / Deformation / Breakage

Unique “Variable Heat Treatment” technology offers superior resistance to torsional stress and breakage compared to other glide path files on the market.



Overcome Calcifications & Difficult Anatomy

A triangular cross-section and rotary motion are highly efficient at overcoming canal impediments.



Simplified Procedure / Reduced Chair Time

Compared to traditional hand filing, Traverse requires fewer instruments for the creation of glide path, leading to easier procedures.



Maintain Original Anatomy / Minimally Invasive

Preserves more natural tooth structure compared to continuously and progressively tapered files. Minimizes unnecessary dentin removal in coronal 1/3.



Superior Rotary Glide Path

Compared to traditional hand filing, Traverse requires fewer instruments leading to easier procedures. Enhance safety by reducing workload and avoiding blockages.



Non-Cutting Tip With Variable Helical Angle Design

Minimizes the chance for ledging or perforation of the canal. Optimized torsional strength makes files less likely to separate and gives clinician greater control.

Ordering Information

818-2146	Traverse Orifice and Glide Path File .06/.13/21mm
818-2147	Traverse Orifice and Glide Path File .06/.13/25mm
818-2148	Traverse Orifice and Glide Path File .06/.13/31mm
818-2196	Traverse Orifice and Glide Path File .06/.18/21mm
818-2197	Traverse Orifice and Glide Path File .06/.18/25mm
818-2198	Traverse Orifice and Glide Path File .06/.18/31mm
818-8256	Traverse Orifice Opener .08/.25/17mm Traverse
818-2156	Traverse Rotary Glide Path File .06/.13/21mm
818-2157	Traverse Rotary Glide Path File .06/.13/25mm
818-2158	Traverse Rotary Glide Path File .06/.13/31mm
818-2186	Traverse Rotary Glide Path File .06/.18/21mm
818-2187	Traverse Rotary Glide Path File .06/.18/25mm
818-2188	Traverse Rotary Glide Path File .06/.18/31mm

