














3M™ Filtek™ Dental Restoratives

A full line of restorative solutions to fit your unique needs.




3M™ Filtek™ Dental Restoratives have been used in more than 2 billion restorations worldwide. Our family of Filtek dental restoratives lets you choose just the right option to fit any clinical need – from highly aesthetic, to highly simple for an easier day.

Properties	Aesthetics & Versatility		Aesthetics & Simplicity		Efficiency & Simplicity	
	3M™ Filtek™ Supreme XTE Universal Restorative	3M™ Filtek™ Supreme Flowable Restorative	3M™ Filtek™ Universal Restorative	3M™ Filtek™ Easy Match Universal Restorative	3M™ Filtek™ One Bulk Fill Restorative	3M™ Filtek™ Bulk Fill Flowable Restorative
 Primary benefits	<ul style="list-style-type: none"> Beautiful aesthetics with natural-looking fluorescence Outstanding strength for anterior and posterior 	<ul style="list-style-type: none"> Virtually no bubbles or “run-on” Ergonomic design Selection of cannula tip sizes (20G, 27G) 	<ul style="list-style-type: none"> Universal opacity for consistent aesthetics Pink Opaquer to mask dark areas Low shrinkage stress 	<ul style="list-style-type: none"> Quick, easy and intuitive shade selection without sacrificing opacity Simplified single-layer restorations, with no blocker needed Naturally-adaptive opacity 	<ul style="list-style-type: none"> Increased opacity for improved aesthetics One-step placement up to 5mm 	<ul style="list-style-type: none"> Virtually no bubbles or “run-on” One-step placement up to 4mm Ergonomic design
 Number of shades	36 (Dentine, Body, Enamel, Translucent)	12	8 (plus Extra White and Pink Opaquer)	3	5	4
 TRUE nanotechnology	Yes	Yes	Yes	Yes	Yes	No
 Handling	Firm	Flows on demand	Creamy	Firm	Sculptable	Adapts easily
 Polish retention	Brilliant					N/A*
 Wear resistance	Excellent					N/A*
 Safe for composite warming and injection molding^{1,2}	Capsules	Syringes	Capsules	Capsules	Capsules	Syringes
						

1. Warm capsules to up to 70°C for up to 1 hour. 2. Warm flowable syringes to up to 70°C for up to 1 hour, up to 25 times. 3. Internal data.

*Filtek Bulk Fill Flowable’s intended use requires allowing at least 2 mm on the occlusal cavosurface for a universal or posterior composite.

We have you covered for every class of direct restoration.

	Aesthetics & Versatility		Aesthetics & Simplicity		Efficiency & Simplicity	
	3M™ Filtek™ Supreme XTE Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ Supreme Flowable Restorative (Incremental up to 2 mm)	3M™ Filtek™ Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ Easy Match Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ One Bulk Fill Restorative (Bulk up to 5 mm)	3M™ Filtek™ Bulk Fill Flowable Restorative (Incremental up to 4 mm)
Class I 	✓	✓	✓	✓	✓ (up to 4 mm)	✓ (Liner/Base)
Class II 	✓	✓	✓	✓	✓ (up to 5 mm)	✓ (Liner/Base)
Class III 	✓	✓	✓	✓		✓
Class IV 	✓	✓	✓	✓		
Class V 	✓	✓	✓	✓		✓
Minimally Invasive	✓	✓	✓	✓	✓	✓
Stress Bearing	✓	✓	✓	✓	✓	

✓ = Recommended use.

www.3M.com.au/Filtek

Request
a demo

