

ULTICLEAN UC600L

Scaler and Air Polisher



Power Meets Comfort



REF: B2-BA260100

Lightweight, anti-slip design

provides a comfortable grip, providing stability and ease of use

Integrated LED light

improves visibility and precision

LED-lit powder container

helps distinguish between powders, streamlining clinical efficiencies

Saves time

with one-touch auto-cleaning (20s and 120s)

Adaptive interface

supports 10 languages, adjustable brightness, and maintenance alerts for optimal performance

Durable handpiece system

offers superior biocompatibility, corrosion, and enhanced wear-resistance

Dual heating system

delivers warm water for both ultrasonic and airstream functions (Ultrasonic scaling temperature: 24°C, 28°C, 32°C, 36°C. Air polishing temperature: 25°C, 30°C, 35°C, 40°C)

Readily accessible, mobile cart available

neatly stores external tubing and cables for a clutter-free workspace

New improved design

provides a quieter airstream (75.8 dB; 9 dB lower than commonly used settings) for an improved patient experience

6 ULTRASONIC MODES FOR CLINICAL VERSATILITY

The Ulticlean UC600L offers robust performance and a wide range of operating modes to meet various clinical needs.

G Mode



Ultra scaling

For rapid and efficient removal of large calculus, reducing procedure time



Soft scaling

For gentle removal of stains and soft deposits, minimising patient discomfort

P Mode



Perio

For effective biofilm removal, protecting soft tissue and minimising trauma



Implant maintenance

Six adjustable levels for precise and safe implant surface cleaning

E Mode



Endo

Precise ultrasonic power optimises root canal preparation and enhances treatment success



Irrigation

Six safety levels ensure effective debris removal while protecting dentin

Order today at
HenrySchein.com.au

 **HENRY SCHEIN®**