

EdgeEndo EdgeFlow Gutta Percha Injecting Needles Gauge and Lengths



Model	Product Code	Gauge	Length
E20G-S	EN-EE20GSIN	20Ga	22mm
E20G-M	EN-EE20GMIN	20Ga	24mm
E20G-L	EN-EE20GLIN	20Ga	28mm
E23G-M	EN-EE23GMIN	23Ga	24mm
E23G-L	EN-EE23GLIN	23Ga	28mm
E25G-M	EN-EE25GMIN	25Ga	24mm
E25G-L	EN-EE25GLIN	25Ga	28mm