


# Testing in your hands with confidence and peace of mind

*Protecting the people you care about  
with SARS-CoV-2 Antigen Self Test Nasal*




Quick and convenient self testing of COVID-19  
everywhere at anytime


**Four simple steps** and the result will be ready after **15 minutes**



Collect nasal sample




Prepare sample



Perform test



Read Result



**Easy-to-use**  
Simple nasal sample collection




**Quick**  
Accurate result in just 15 minutes



**Global leader in diagnostics**  
125 year history, operating in 150 countries



# The simple, convenient way to manage health credentials with NAVIFY® Pass



**NAVIFY® Pass is a free digital COVID-19 solution** to allow you to record and share your diagnostic test results. This can be used for your personal record of self testing using the Roche SARS-CoV-2 Antigen Self Test Nasal.

After performing the Self Test, simply upload the test result to NAVIFY® Pass and easily display your COVID-19 health status to others where needed.

## The NAVIFY® Pass is simple to use:

- 1 Download NAVIFY® Pass from app store
- 2 Register NAVIFY® Pass account in the mobile app
- 3 Perform Self Test
- 4 Scan data matrix on test device
- 5 Fill out questionnaire and submit test result
- 6 Share COVID-19 health-data where needed

**NAVIFY® Pass is free** and can be downloaded from the App Store and Google Play.



**Technical Support Call Hours: 9am-7pm AEST / 9am-8pm AEDT, 7 days per week 1800 497 069**

NAVIFY is trademark of Roche. All other product names and trademarks are the property of their respective owners.  
 © 2021 Roche Diagnostics. Roche Diagnostics Australia Pty Limited. 2 Julius Avenue North Ryde NSW 2113 Australia.  
 ABN 29 003 001 205. Phone: 02 9860 2222 E-mail: australia.dia\_cardiology@roche.com.au MC-AU-01413

**NAVIFY®**