

# HENRY SCHEIN NEEDLES





## HENRY SCHEIN NEEDLES



CODE	DESCRIPTION
HS-5723808	<a href="#">HS Standard Needle 27G Short Orange 25mm Box of 100</a>
HS-5723809	<a href="#">HS Standard Needle 30G Short Blue 25mm Box of 100</a>
HS-5723816	<a href="#">HS Standard Needle 27G X-Long Yellow 40mm Box of 100</a>
HS-5723807	<a href="#">HS Standard Needle 27G Long Yellow 35mm Box of 100</a>
HS-5723811	<a href="#">HS Standard Needle 30G X-Short Purple 12mm Box of 100</a>

Sterile, offers superior precision, comfort, and reliability for all your clinical needs.

- **Plastic Hub with Indicating Mark:** Easily identify bevel position thanks to the clear marking on the hub, ensuring accurate and precise needle placement during procedures.
- **Round Cutting Edge:** Our round cutting edge design protects the patient's gum, minimizing discomfort while maintaining the integrity of the tissue.
- **Double Silicone Coating:** Each needle features a double silicone coating that reduces tissue trauma, making injections smoother and more comfortable for the patient.
- **Colour-Coded Seals & Hubs:** For quick and easy identification, seals and hubs are color-coded according to diameter, ensuring you always grab the correct size for each procedure.



[Contact us](#)

 HENRY SCHEIN®