

PROPHYflex™ 4
Nex-gen air polisher.



Two handpieces – one streamlined device.

Ready for a daily prophylaxis treatments that's easier, quicker, more thorough, and even enjoyable? The new KaVo PROPHYflex 4 air polisher checks all the boxes. It's feather-light, perfectly balanced, and ergonomically designed for a more comfortable experience.

Exceptional cleaning. Enhanced visibility.

Less powder mist thanks to the highly precise powder jet. Plus, the narrow shape of the cannula provides improved visibility during supragingival cleaning.

360 of comfort.

Work comfortably and access hard-to-reach spots thanks to the 360° rotation grip sleeve.

Adjusted to fit you.

The ultra-light, perfectly balanced handpiece has individually adjustable lengths for comfortable handling.



Comfort and Flexibility.

"The new PROPHYflex 4 has excelled in my daily work – it fits perfectly in the hand and can be individually configured through the two grip sleeve lengths. The spray mist can also be easily controlled, which is something my patients like as well."

Tanja Wicharz

Prophylaxis expert and practice instructor

Because every hand, every treatment, and every patient is different, KaVo PROPHYflex 4 can be individually adjusted to each user. Whether you're working in subgingival or supragingival regions, you get to decide the amount of powder used to eliminate discoloration and how long the handpiece should be.

More precision for pocket areas.

The optional Perio Kit provides a focused subgingival treatment, gently cleans bridges, implants and crowns.

Tailored to you.

Two grip sleeves with different lengths to fit your hand perfectly. Also ideal when several practitioners work with the same handpiece.

Simplicity and safety for success.

The powder container can be filled in advance and secured to simplify your workflow and avoid interruptions during treatment.

Three levels of adjustment.

Level 1 delivers powder-less cleaning of the teeth to remove any powder residue. Levels 2 and 3 provide individually adjustable cleaning performances depending on the patient and indication.


Flexible.

Use with KaVo MULTIflex coupler for the ultimate in flexibility. Easily move from room to room.



Before, during, and after treatment.

For the best possible infection control, KaVo PROPHYflex 4 is seamlessly designed so you can easily disassemble, clean, disinfect, and sterilize the device each time.




Less powder mist.



You will be delighted with the adjustable powder quantities and constant, focused powder jet which brings precise cleaning performance and minimized powder mist formation, for an optimal hygiene routine.

Tough on stains. Gentle on teeth.



Energy is evenly distributed along the jet surface for more effective cleaning performance and greater protection of the tooth structure.

* PROPHYpearls abrasion compared to a competitor's air polishing device using calcium carbonate. Shown on KaVo Kerr SonicFill™ composite.



Control with adjustable power settings.

With three adjustment settings, you have the power to decide the cleaning performance, from powder-free to the maximum cleaning mode.



New color options.



PROPHYflex 4 doesn't just match the shape of your hand, it also matches your style. Choose from a Flamingo pink, Wave blue, or Lime green color options.

Save time with sealable containers.



Seal-able powder containers are ideal for preparing for treatments where the choice of powder needs to match the indication. Simply attach the pre-filled container and go.

PROPHYflex 4 cleaning powders.


PROPHYpearls™, PROPHYflex™ Powder, and PROPHYflex™ Perio Powder are designed for the KaVo PROPHYflex 4 and allow for a thorough and gentle teeth cleaning.



KaVo PROPHYpearls:

For efficient clean and gentle prophylaxis.


- **Indication:** effective elimination of discolouration, plaque and biofilm.
- **Active ingredient:** calcium carbonate at a particle size of \varnothing 60-70 μ m, biodegradable, no salt, sugar or gluten.
- **Highlights:** round shape and a uniform size for optimal cleaning with minimal abrasion, four flavors and two neutral variants.



KaVo PROPHYflex Powder:

For speedy cleaning.

- **Indication:** highly effective elimination of discoloration, plaque and biofilm.
- **Active ingredient:** sodium bicarbonate at a particle size of \varnothing 60-70 μ m, water-soluble.
- **Highlights:** simple, fast and effective treatment with four flavors.



Perio Powder for subgingival procedures.

With the PROPHYflex 4, you no longer need to change the device when switching from supragingival to subgingival treatments or implant cleaning. Simply attach the perio tip, fill with Perio Powder, and choose your desired performance level. It's that simple.


- Safe, targeted pocket cleaning thanks to the stable, thin perio tip with depth marking
- Excellent accessibility thanks to good visibility and a slim shape.
- Hygienic and multi-application: each tip can be sterilized up to 10 cycles*

*Data on file



KaVo PROPHYflex Perio Powder: for subgingival and supragingival treatment.

- **Indication:** efficient biofilm management in periodontal pockets and cleaning of crowns, bridges and implants.
- **Active ingredient:** fine glycine powder (18-22 μm), protects gums.
- **Highlights:** gentle and not too abrasive, suitable for efficient cleaning of pockets as well as supragingival applications on implants.





Handpieces

KaVo has always been the leader in creating innovative solutions for dental practitioners. Our vast line of quality handpieces showcase our attention to your level of care while delivering performance that lasts.



Treatment Units

Beautiful lines, patient comfort and simple operation are just a few of the benefits to the line of KaVo treatment units. Everything you need to perform any procedure –all in one solution.



Imaging Solutions

Designed with ease-of-use for all clinicians in mind, KaVo now offers dependable and consistent imaging solutions that provide vital information to support accurate diagnosis and predictable treatment planning.

The products, equipment and services illustrated and described in this brochure reflects knowledge at the time of printing. KaVo Dental accepts no liability for any deviation from the illustrations in terms of color or shape, or any errors or print errors, and retains the right to make changes to the brochure at any time. Full or partial reprinting is only permitted with permission from KaVo Dental. PROPHYflex, MULTIflex, PROPHYpowder, PROPHYpearls, are registered trademarks of KaVo Dental in the United States and/or other countries.

© 2019 KaVo Dental. All rights reserved. KV00576/A8.19

KaVo | 11727 Fruehauf Drive | Charlotte NC 28273 | USA | 1-888-ASK-KAVO

kavo.com