

Komet Dental Clinical Burs



Tungsten Carbide Instruments

	H1 204 005-027 314 005-027 316 006-029
	H1SEM 204 010, 012, 014, 016, 018-031 205 010, 012, 014, 016, 018-027
	H1SOS 320 023
	H33LSOS 320 016
	H7 204 008, 010, 014 313 008 314 006-018
	H21 204 008-014 314 008-014
	H21R 314 010, 012, 014
	H23R 204 012, 016 314 010, 012, 016
	H31 204 010, 012, 014, 016 314 008-018 316 010, 012, 014
	H33 104 009, 010, 012, 016, 021 314 009, 010, 012, 016 316 010, 012, 016, 021
	H33L 104 009, 010, 012, 016, 021 106 016 (70mm HP Length) 314 009, 010, 012, 016, 021 316 010, 012, 016
	H33R 316 016
	H141 104 010-050 106 023, 031 (70mm HP Length) 205 010-040
	H254E 104 012 204 012 314 012
	H161 104 016 314 016
	H162 104 016 (H162SXL option available in XL Friction Grip, 31mm Length) 204 016 205 016 314 016
	H379 204 014, 018, 023 314 014, 018, 023
	H4MC (& L/XL/XXL) 314 010, 012 (also L) 314 014 XL/XXL

Tungsten Carbide Q-Finisher Instruments

	H135Q 314 014 Palatal/Occlusal (Non-Cutting Tip)
	H48LQ 314 012 Labial (Non-Cutting Tip)
	H246Q 314 009 Occlusal (With Cutting Tip)

Ceramic Instruments

	K1SM 204 008-016, 018, 021, 023, 027 205 010-016, 018, 021, 023
	K25 1ACR 104 060 Also available in Tungsten option

Diamond Instruments

All shank/size options indicated below are for Medium grits.

	801 ○ ● ● ● ● ● 313 010, 014, 018 314 006-023, 029, 035 316 014, 023, 029 204 009-023, 029, 035
	835KR ● ● ● ● ● ● 314 007-016
	368 ○ ● ● ● ● ● 204 023 314 016, 021, 023
	379 ○ ● ● ● ● ● 314 014, 018, 023
	836KR ● ● ● ● ● ● 314 010-018
	837KR ● ● ● ● ● ● 314 010, 012, 014
	838 ● ● ● ● ● ● 314 008, 009, 010, 012
	847KR ○ ● ● ● ● ● ● 314 014, 016, 023
	848 ● ● ● ● ● ● 314 014, 016, 018, 023
	850 ● ● ● ● ● ● 314 012, 014, 016, 018, 023
	858 ○ ● ● ● ● ● ● 104 014 314 010, 014, 016
	830 ● ● ● ● ● ● 313 012 314 007-010, 012, 014, 016
	862 ○ ● ● ● ● ● ● 204 014, 018 314 010, 012, 014, 016, 021
	878 ● ● ● ● ● ● 313 012 314 010, 012, 014, 016
	852 ○ ● ● ● ● ● ● 314 012, 037
	881 ● ● ● ● ● ● ● 314 010, 012, 014, 016
	856 ● ● ● ● ● ● ● 314 012, 014, 016, 018, 021
	806 ● ● ● ● ● ● ● 314 009-018
	888 ● ● ● ● ● ● ● 313 012 314 012
	909 ● ● ● ● ● ● ● 204 040, 055, 065 314 040, 055
	6844 ● ● ● ● ● ● ● 314 014, 016 2-Gritted Veneer Bur
	8392 ● ● ● ● ● ● ● 314 016 Interproximal Trimmer
	370 ● ● ● ● ● ● ● 314 030, 035 OccluShaper
	DM ● ● ● ● ● ● ● 314 009, Available in 4 depths: 0.5, 1.0, 1.5, 2.0mm PrepMarker
	4ZR ● ● ● ● ● ● ● 314 012, 014 Zirconia Crown Cutter

Diamond Strips

All options are 148mm long

Medium Grit = 45µm
Fine Grit = 30µm
Ultra Fine Grit = 15µm

	DS25/37 (& F/EF)* ● ● ● ● ● ● ●
	WS25/37 (& F/EF)* ● ● ● ● ● ● ●

* Assorted "A" grit packs available in both DS25 & WS25 (2.5mm) & DS37 & WS37 (3.7mm) width options. New DS60 (6.0mm) wide option available in Medium & Fine grits only.

Endodontic Instruments

	H269QGK 314 012, 016 315 016 (21mm Length) EndoGuard (Endo-Z)
	H1SML31/34 205 004, 006, 008, 010, 012, 014 (in 31 or 34mm Length Options) EndoTracer
	OP08L (Taper 08) 204 025 (19mm Length Only) Opener
	OP10L (Taper 10) 204 030 (15mm or 19mm Lengths)
	17521/25 (Hand File, Taper 02) 21mm Length in: 006, 008, 010, 012, 015 25mm Length in: 006, 008, 010, 012, 015 Patency (Heat-treated K-File)
	PG03L (Taper 03) 204 015, 020 (in 21, 25, 31mm Length Options) PathGlider
	R06L (Taper 06) 204 025, 030, 035, 040, 045, 050 (in 21, 25, 31mm Length Options) Reziflow Reciprocating NiTi File

Steel Instruments

	1 204 005-031 205 010-023
	41 204 010-023
	191 204 090, 100, 120, 140, 160, 180 Endo Pulp Bur
	9816 Picture features exaggerated cutting edge Interproximal Separating Strip

Information

Diamond grit sizes

○ ultra-fine 8µm	● medium 107µm*
● extra-fine 25µm	● coarse 151µm*
● fine 46µm	● super-coarse 181µm*
	● extra-coarse 230µm

* Depending on bur shape/size, grit sizes may vary slightly from specified values.

Shank Types

	314 Friction Grip (FG) Ø 1.60mm
	313 Friction Grip Short-Paedo (FGS) Ø 1.60mm
	316 Friction Grip Extra Long (FGXL) Ø 1.60mm
	204 Right Angle (RA) Ø 2.35mm
	205 Right Angle Long (RAL) Ø 2.35mm
	104 Handpiece (HP) Ø 2.35mm
	320 Surgical Length Handpiece (HP) with Notched Neck Ø 2.35mm