


NeoLucent Plus
Ceramic Bracket System



NeoLucent Plus Ceramic Bracket System

The *NeoLucent® Plus™* Low-Profile Ceramic Bracket System delivers superb aesthetics by blending with the natural tooth tone for a virtually invisible appearance. With improved performance due to a *Crunch COAT™* Base, more anatomically contoured base radii, and improved materials for a smoother finish and color consistency, *NeoLucent Plus* Ceramic Brackets provide the definitive performance you demand.

Henry Schein® Orthodontics™ is dedicated to providing world class products and solutions to meet your every orthodontic need. As a full-service partner, we offer the tools to help expand your practice and diversify treatment options. So when your patients desire a cosmetic alternative, we have the right choice for you, and them.

PATIENT SATISFACTION

Invisible Appearance

Crafted to the highest aesthetic standards

Translucent, natural appearance

from 99.9% pure Polycrystalline Alumina

Stain resistant

throughout treatment

COMFORT

Designed with your patient's comfort in mind

Low profile

Smooth, polished surface

Rounded Corners

WHAT'S IN THE PLUS?

Crunch Coat Base


Base surface is proprietary to *Henry Schein Orthodontics*

Bond Strength Comparable to Mesh Pad

Fine alumina particles bonded to the bracket base provide a large surface area for adhesive bonding

Consistent Debonding

Because debonding occurs near the *Crunch COAT* adhesive interface, there is minimal distress to the enamel surface



CLINICAL PERFORMANCE

Accurate placement achieved with an easily removable color-coded ID mark

Strong bonding and predictable debonding with a *Crunch COAT Base*

Superior Sliding Mechanics via a precision-cut archwire slot


Non-allergenic for sensitive patients as brackets do not contain nickel, chromium or other metals

Easy Ligation with ample tie wing clearance

Routine Protocol with no special instruments or adhesives required


SINGLE HOOK

More Comfort for the Patient, less bulky than universal hooks




COMPATIBLE IN/OUT WITH MAESTRO LOW-PROFILE BRACKET SYSTEM AND NEOCRYSTAL SAPPHIRE BRACKET SYSTEM

For aesthetic anterior and metal posterior bracketing



Maestro® Bracket



NeoLucent Plus Bracket



NeoCrystal™ Bracket

The mission of *Henry Schein Orthodontics* is to provide state-of-the-art orthodontic products and innovative clinical solutions that enable our customers to offer exceptional patient care, while expanding the scope and profitability of their practices.



P: 1300 ORDERZ (673 379)
E: sales@naol.com.au

