

# The Original CCLAD System Technology

Since 1997

*"The Wand® is one of those products that I can say honestly changed my life, and one I can't practice without."*

Michael DiTolla, DDS, FAGD

*"One of the first things our patients experience is the computerized injection using the Wand®. Compliments follow, as many have never had such a wonderful experience with a dental injection."*

Alan N. Rembos, DDS

*"The Wand® increases patient loyalty, drives referrals, and helps define my brand. It's one of the top 10 reasons for my success, and is an instrument I can't practice without."*

Cynthia K. Brattesani, DDS

*"The Wand® is better than the traditional syringe for patient comfort and it causes much less anxiety. The Wand differentiates our practice. Patients love high-tech dentistry and it gives them a much better overall experience."*


David Hornbrook, DDS, FAGD

# the Wand®

The Original CCLAD System Technology

Since 1997

MORE THAN  
3,000,000  
HANDPIECES SOLD PER YEAR  
WORLDWIDE



## A New Marketing Opportunity

DIFFERENTIATES THE PRACTICE • INCREASES REFERRALS • DRIVES LOYALTY

## THE GLOBAL LEADER IN CCLAD System Technology

COMPUTER CONTROLLED LOCAL ANESTHETIC DELIVERY

# the Wand®

Computer-Assisted Anesthesia System

**thewand.com • 800-862-1125**

Wand Dental, Inc. • 425 Eagle Rock Avenue • Roseland, NJ 07068



# Dentist-Controlled Computer-Assisted

## COMPUTER-REGULATED FLOW RATES

- All injection techniques can be performed **below the patient's pain threshold**
- Patented technology **precisely controls flow rate and pressure of anesthesia**
- 3 Speeds are **tailored for each injection type**

## DYNAMIC PRESSURE SENSING DPS®

- Enables successful, virtually painless **single-tooth anesthesia with no collateral numbness**
- Patented technology monitors exit pressure to **ensure optimal needle position is maintained**
- Visual and audible feedback **aids accurate needle placement** for STA/PDL injection

## THE WAND® HANDPIECE

- Unparalleled tactile **control**
- **Reduces anxiety—looks nothing like a syringe**
- Eliminates needle deflection, **reducing number of missed blocks**
- Small intra-oral profile **increases clinician visibility**
- Easily concealable and flexible size to **improve access**

# the Wand®

Eliminates the Pain &  
Anxiety of Injections



# More Comfort, More Options...

EVERY INJECTION IS BETTER WITH THE WAND®

## STA-INTRALIGAMENTARY (Modified PDL)

- About 80% of dental procedures are performed on a single tooth.
- **Reduces need for mandibular blocks**
- Increases patient comfort
- **Eliminates collateral numbness**
- Immediate onset of anesthesia



## SUPRA-PERIOSTEAL (Infiltrations)

- Controlled flow rate ensures the fluid pressure of the injection is below the patient's pain threshold



## INFERIOR ALVEOLAR (Blocks)

- The Wand® is held like a pen, enabling bi-rotational techniques—eliminating needle deflection, reducing missed blocks
- Can provide multi-cartridge anesthesia with a single injection



## PALATALS (P-ASA & AMSA)

- Patented flow rate and pressure settings allow for profound pulpal anesthesia while greatly improving patient comfort
- P-ASA & AMSA techniques allow for immediate smile line assessment (no collateral anesthesia)

