



Screw guide





EXTERNAL HEX UNIVERSAL

compatible with
Nobel Brånemark®

SOCKET
UniGrip®



DESS

> DENTAL SMART SOLUTIONS

	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.027	Screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.028	Screw UG RP	2,53	0°	7,35	1,42	3,12	2,81	M2,0	35 Ncm
	19.029	Screw UG RP	2,85	0°	7,10	1,40	2,70	3,00	M2,5	35 Ncm
	19.227	DLC screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.228	DLC screw UG RP	2,53	0°	7,35	1,42	3,12	2,81	M2,0	35 Ncm
	19.229	DLC screw UG RP	2,85	0°	7,10	1,40	2,70	3,00	M2,5	35 Ncm
	19.627	TIN screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.628	TIN screw UG RP	2,53	0°	7,35	1,42	3,12	2,81	M2,0	35 Ncm
	19.629	TIN screw UG RP	2,85	0°	7,10	1,40	2,70	3,00	M2,5	35 Ncm




TRI-LOBE

compatible with
NobelReplace Select™

SOCKET
UniGrip®



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.004	Screw UG NP	2,45	0°	8,30	1,70	2,60	2,75	M1,6	35 Ncm
	19.005	Screw UG RP & WP	2,45	0°	9,90	1,50	3,12	2,81	M2,0	35 Ncm
	19.204	DLC screw UG NP	2,45	0°	8,30	1,70	2,70	3,00	M2,5	35 Ncm
	19.205	DLC screw UG RP & WP	2,45	0°	9,90	1,50	2,60	2,75	M1,6	35 Ncm
	19.604	TIN screw UG NP	2,45	0°	8,30	1,70	3,12	2,81	M2,0	35 Ncm
	19.605	TIN screw UG RP & WP	2,45	0°	9,90	1,50	2,70	3,00	M2,5	35 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



ACTIVE HEX

compatible with
NobelActive® / Replace® CC

SOCKET
UniGrip®



DESS

> DENTAL SMART SOLUTIONS

	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.027	Screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.030	Screw UG RP	2,53	54°	7,35	1,57	2,97	2,81	M2,0	35 Ncm
	19.068	Screw UG 3.0	1,82	55°	6,72	1,82	2,45	2,45	M2,5	15 Ncm
	19.227	DLC screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.230	DLC screw UG RP	2,53	54°	7,35	1,57	2,97	2,81	M2,0	35 Ncm
	19.268	DLC screw UG 3.0	1,82	55°	6,72	1,82	2,45	2,45	M2,5	15 Ncm
	19.627	TIN screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.630	TIN screw UG RP	2,53	54°	7,35	1,57	2,97	2,81	M2,0	35 Ncm
	19.668	TIN screw UG 3.0	1,82	55°	6,72	1,82	2,45	2,45	M2,5	15 Ncm




OCTAGON

compatible with
Straumann® Tissue level & synOcta®

SOCKET
Torx® 06



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.009	Screw Torx® on abutment	2,50	45°	4,50	1,20	0	3,30	M2,0	15 Ncm
	19.026	Screw Torx®	2,58	75°	6,70	3,20	1,68	2,00	M2,0	35 Ncm
	19.070	Screw Torx® NNC	2,10	75°	7,10	1,85	3,35	1,70	M1,6	35 Ncm
	19.609	TIN screw Torx® on abutment	2,50	45°	4,50	1,20	0	3,30	M2,0	15 Ncm
	19.626	TIN screw Torx®	2,58	75°	6,70	3,20	1,68	2,00	M2,0	35 Ncm
	19.670	TIN screw Torx® NNC	2,10	75°	7,10	1,85	3,35	1,70	M1,6	35 Ncm



A. HEAD DIAMETER (mm)

B. HEAD BASE ANGLE

C. OVERALL LENGTH (mm)

D. HEAD LENGTH (mm)

E. DISTANCE FROM
HEAD TO THREAD (mm)

F. THREAD LENGTH (mm)

G. THREAD KIND



CONICAL BL

compatible with
Straumann® Bone level

SOCKET
Torx® 06



DESS

> DENTAL SMART SOLUTIONS

	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.032	Screw Torx® NP	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.033	Screw Torx® RP	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.037	Screw Torx® NC/RC	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.098	Screw for SRA	2,20	55°	3,65	1,65	0,20	1,70	M1,4	15 Ncm
	19.632	TIN screw Torx® NP	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.633	TIN screw Torx® RP	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.637	TIN screw Torx® NC/RC	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.698	TIN screw for SRA	2,20	55°	3,65	1,65	0,20	1,70	M1,4	15 Ncm



CONICAL BLX

compatible with
Straumann® BLX

SOCKET
Torx® 06



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.054	Screw Torx®	2,10	90°	6,10	1,02	4,26	1,84	M1,6	35 Ncm
	19.700	TIN screw Torx®	2,10	90°	6,10	1,02	4,26	1,84	M1,6	35 Ncm




NEO GM

compatible with
Neodent® Grand Morse

SOCKET
GM



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.051	Screw UG®	1,86	60°	9,75	1,50	7,65	2,15	M1,6	20 Ncm



A. HEAD DIAMETER (mm)
 B. HEAD BASE ANGLE
 C. OVERALL LENGTH (mm)
 D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
 F. THREAD LENGTH (mm)
 G. THREAD KIND



INTERNAL HEX CONIC

compatible with
Astra Tech Osseospeed™

SOCKET
HEX 1,27mm



DESS

> DENTAL SMART SOLUTIONS

	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.015	Screw hex NP (3.0)	2,10	30°	7,50	2,00	3,00	2,50	M1,4	15 Ncm
	19.016	Screw hex RP	2,35	0°	8,25	2,00	3,22	3,02	M1,6	20 Ncm
	19.017	Screw hex WP	2,35	0°	8,30	2,00	3,80	3,50	M2,0	25 Ncm
	19.018	Screw on abutment	2,30	0°	4,80	1,90	0,00	2,90	M1,4	15 Ncm
	19.615	TIN screw hex NP (3.0)	2,10	30°	7,50	2,00	3,00	2,50	M1,4	15 Ncm
	19.616	TIN screw hex RP	2,35	0°	8,25	2,00	3,22	3,02	M1,6	20 Ncm
	19.617	TIN screw hex WP	2,35	0°	8,30	2,00	4,80	3,50	M2,0	25 Ncm
	19.618	TIN screw on abutment hex	2,30	0°	4,80	1,90	0,00	2,90	M1,4	15 Ncm




CONIC EVO

compatible with
Astra Tech implant system™ EV

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.058	Screw hex on Uniabutment EV	2,35	18°	4,49	1,71	0,00	2,50	M1,8	15 Ncm
	19.059	Screw hex EV 3.6	2,15	60°	7,46	2,12	1,34	4,00	M1,6	25 Ncm
	19.060	Screw hex EV 4.2	2,35	40°	7,64	1,85	1,79	4,00	M1,8	25 Ncm
	19.061	Screw hex EV 4.8/5.4	2,35	0°	7,60	1,60	2,00	4,00	M2,0	25 Ncm
	19.071	Screw hex EV 3.0	2,15	52°	7,68	1,60	3,83	3,85	M1,4	25 Ncm
	19.658	TIN screw hex on Uniabutment EV	2,35	18°	4,49	1,71	0,00	2,50	M1,8	15 Ncm
	19.659	TIN screw hex EV 3.6	2,15	60°	7,46	2,12	1,34	4,00	M1,6	25 Ncm
	19.660	TIN screw hex EV 4.2	2,35	40°	7,64	1,85	1,79	4,00	M1,8	25 Ncm
	19.661	Screw hex EV 4.8/5.4	2,35	0°	7,60	1,60	2,00	4,00	M2,0	25 Ncm



A. HEAD DIAMETER (mm)

B. HEAD BASE ANGLE

C. OVERALL LENGTH (mm)

D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)

F. THREAD LENGTH (mm)

G. THREAD KIND



INTERNAL HEX FD

compatible with
Friadent Xive®

SOCKET
HEX 1,22mm



DESS

> DENTAL SMART SOLUTIONS



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.024	Screw hex	2,10	45°	9,80	3,06	3,29	3,45	M1,6	25 Ncm



19.624	TIN screw hex	2,10	45°	9,80	3,06	3,29	3,45	M1,6	25 Ncm
--------	---------------	------	-----	------	------	------	------	------	--------



INTERNAL ANK

compatible with
Dentsply Ankylos® C/X

SOCKET
HEX 1,00mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.048	Screw hex	2,35	60°	7,87	2,37	3,40	3,10	M1,8	15 Ncm



19.050	Screw hex	2,35	60°	9,37	2,02	3,35	3,40	M1,8	15 Ncm
--------	-----------	------	-----	------	------	------	------	------	--------



19.648	TIN screw hex	2,35	60°	7,87	2,37	2,40	3,10	M1,8	15 Ncm
--------	---------------	------	-----	------	------	------	------	------	--------



EXTERNAL HEX USA

compatible with
3i® Osseotite®

SOCKET
HEX 1,20mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.034	Screw hex	2,55	0°	7,85	1,65	0,70	5,50	M2,0	20 Ncm



19.634	TIN screw hex	2,55	0°	7,85	1,65	0,70	5,50	M2,0	20 Ncm
--------	---------------	------	----	------	------	------	------	------	--------



INTERNAL HEX "CLICK"

compatible with
3i Certain®


SOCKET
HEX 1,20mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.011	Screw hex	2,00	0°	8,40	1,90	3,80	2,70	M1,6	20 Ncm



19.611	TIN screw hex	2,00	0°	8,40	1,90	3,80	2,70	M1,6	20 Ncm
--------	---------------	------	----	------	------	------	------	------	--------



A. HEAD DIAMETER (mm)

B. HEAD BASE ANGLE

C. OVERALL LENGTH (mm)

D. HEAD LENGTH (mm)

E. DISTANCE FROM
HEAD TO THREAD (mm)

F. THREAD LENGTH (mm)

G. THREAD KIND



INTERNAL HEX USA

compatible with
Zimmer Screw-Vent®

SOCKET
HEX 1,27mm



DESS

> DENTAL SMART SOLUTIONS



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.012	Screw hex	2,30	19°	8,20	2,55	2,98	2,60	UNF 1-72	30 Ncm
19.042	Screw hex	1,95	19°	8,20	2,49	2,95	2,76	UNF 1-72	30 Ncm
19.642	TIN screw hex	1,95	19°	8,20	2,49	2,95	2,76	UNF 1-72	30 Ncm



INTERNAL HEX MI

compatible with
MIS® Seven

SOCKET
HEX 1,27mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.056	Screw hex 3.0	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 Ncm
19.081	Screw hex SP/WP	2,20	28°	7,48	2,41	2,82	2,25	UNF 1-72	30 Ncm
19.656	TIN screw hex 3.0	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 Ncm
19.681	TIN screw hex SP/WP	2,20	28°	7,48	2,41	2,82	2,25	UNF 1-72	30 Ncm
19.110	Screw hex RP	2,00	0°	3,45	1,35	1,35	2,10	M1,4	15 Ncm



CONIC MC

compatible with
MIS® C1 internal

SOCKET
HEX 1,27mm




CONIC MV

compatible with
MIS® V3

SOCKET
HEX 1,27mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.056	Screw hex NP	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 Ncm
19.081	Screw hex SP/WP	2,20	28°	7,48	2,41	2,82	2,25	UNF 1-72	30 Ncm
19.656	TIN screw hex NP	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 Ncm
19.681	TIN screw hex SP/WP	2,20	28°	7,48	2,41	2,82	2,25	UNF 1-72	30 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



EXTERNAL HEX BH

compatible with
BioHorizons® External

SOCKET
HEX 1,27mm



DESS

> DENTAL SMART SOLUTIONS

	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.013	Screw hex NP	2,50	0°	7,00	1,00	2,10	3,90	UNF 1-72	30 Ncm
	19.014	Screw hex RP/WP	2,65	0°	7,10	1,00	2,10	4,00	M2,0	30 Ncm
	19.613	TIN screw hex NP	2,50	0°	7,00	1,00	2,10	3,90	UNF 1-72	30 Ncm
	19.614	TIN screw hex RP/WP	2,65	0°	7,10	1,00	2,10	4,00	M2,0	30 Ncm



INTERNAL HEX BH

compatible with
BioHorizons® internal

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.043	Screw hex	2,40	45°	7,24	2,42	2,36	2,45	UNF 1-72	30 Ncm
	19.643	TIN screw hex	2,40	45°	7,24	2,42	2,36	2,45	UNF 1-72	30 Ncm




INTERNAL CAM

compatible with
Camlog®

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.044	Screw hex NP	2,10	45°	10,45	2,35	4,60	3,50	M1,6	20 Ncm
	19.045	Screw hex WP	2,50	45°	10,45	2,35	4,60	3,50	M2,0	20 Ncm
	19.644	TIN screw hex NP	2,10	45°	10,45	2,35	4,60	3,50	M1,6	20 Ncm
	19.645	TIN screw hex WP	2,50	45°	10,45	2,35	4,60	3,50	M2,0	20 Ncm



A. HEAD DIAMETER (mm)

B. HEAD BASE ANGLE

C. OVERALL LENGTH (mm)

D. HEAD LENGTH (mm)

E. DISTANCE FROM
HEAD TO THREAD (mm)

F. THREAD LENGTH (mm)

G. THREAD KIND



INTERNAL CON SOCKET
HEX 1,27mm
compatible with
Conelog®



DESS
> DENTAL SMART SOLUTIONS

	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.074	Screw hex 3.3/3.8/4.3	2,26	30°	8,87	2,08	3,99	2,80	M1,6	20 Ncm
	19.076	Screw hex 5.0	2,67	30°	8,87	2,08	3,99	2,80	M2,0	20 Ncm
	19.674	TIN screw hex 3.3/3.8/4.3	2,26	30°	8,87	2,08	3,99	2,80	M1,6	20 Ncm
	19.676	TIN screw hex 5.0	2,67	30°	8,87	2,08	3,99	2,80	M2,0	20 Ncm



CONIC ANYR SOCKET
HEX 1,20mm
compatible with
Megagen AnyRidge®




	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.057	Screw hex	2,10	0°	8,00	2,00	2,78	3,22	M1,8	35 Ncm
	19.657	TIN screw hex	2,10	0°	8,00	2,00	2,78	3,22	M1,8	35 Ncm



CONIC AXIO SOCKET
HEX 1,18mm
compatible with
Anthogyr AXIOM®



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.062	Screw hex Multi-Unit	2,18	0°	4,45	2,00	0,00	2,45	M1,4	15 Ncm
	19.077	Screw hex	2,22	70°	7,90	2,40	3,30	2,20	M1,6	25 Ncm
	19.662	TIN screw hex Multi-Unit	2,18	0°	4,45	2,00	0,00	2,45	M1,4	15 Ncm
	19.677	TIN screw hex	2,22	70°	7,90	2,40	3,30	2,20	M1,6	25 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



CONIC BIO

compatible with
Biotech KONTACT®

SOCKET
HEX 1,22mm



DESS

> DENTAL SMART SOLUTIONS

	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.063	Screw hex NrW	2,03	30°	4,50	1,80	0,00	2,70	M1,6	20 Ncm
	19.064	Screw hex Std	1,93	45°	4,50	1,95	0,00	2,55	M1,4	20 Ncm
	19.078	Screw hex	1,93	45°	7,58	1,56	3,27	2,75	M1,4	20 Ncm
	19.663	TIN screw hex NrW	2,03	30°	4,50	1,80	0,00	2,70	M1,6	20 Ncm
	19.664	TIN screw hex Std	1,93	45°	4,50	1,95	0,00	2,55	M1,4	20 Ncm
	19.678	TIN screw	1,93	45°	7,58	3,27	3,27	2,75	M1,4	20 Ncm



CONIC DENT

compatible with
Dentium Superline™ & Implantium®

SOCKET
HEX 1,28mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.066	Screw hex MUA	2,30	45°	5,00	1,93	0,57	2,50	M1,8	10 Ncm
	19.090	Screw hex S-Line	2,30	33°	10,2	2,05	4,85	3,30	M2,0	30 Ncm
	19.666	TIN screw hex MUA	2,30	45°	5,00	1,93	0,57	2,50	M1,8	10 Ncm
	19.690	TIN screw hex S-Line	2,30	33°	10,2	2,05	4,85	3,30	M2,0	30 Ncm




CONIC IC

compatible with
Medentis ICX®

SOCKET
HEX 1,40mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.080	Screw hex	2,47	75°	9,10	2,97	4,13	2,00	M1,6	30 Ncm
	19.680	TIN screw hex	2,47	75°	9,10	2,97	4,13	2,00	M1,6	30 Ncm



A. HEAD DIAMETER (mm)
 B. HEAD BASE ANGLE
 C. OVERALL LENGTH (mm)
 D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
 F. THREAD LENGTH (mm)
 G. THREAD KIND



CONIC INK

compatible with
GlobalID IN-KONE®

SOCKET
HEX 1,22mm



DESS

> DENTAL SMART SOLUTIONS

	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.035	Screw hex Straight	2,22	45°	5,00	2,10	0,00	2,90	M1,8	15 Ncm
	19.036	Screw hex Angled	2,10	45°	4,80	2,10	0,00	2,70	M1,4	15 Ncm
	19.086	Screw hex	2,17	31°	8,45	2,24	4,21	2,00	M1,8	15 Ncm
	19.635	TIN screw hex Straight	2,22	45°	5,00	2,10	0,00	2,90	M1,8	15 Ncm
	19.636	TIN screw hex Angled	2,10	45°	4,80	2,10	0,00	2,70	M1,4	15 Ncm
	19.686	TIN screw	2,17	31°	8,45	2,24	4,21	2,00	M1,8	15 Ncm



CONIC OSS

compatible with
Osstem® TS / HIOSSEN® ET

SOCKET
HEX 1,22mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.072	Screw hex Mini	2,22	45°	18,18	2,67	5,41	2,10	M1,6	20 Ncm
	19.073	Screw hex Standard	2,33	45°	8,35	2,60	2,95	2,80	M2,0	30 Ncm
	19.272	TIN screw hex Mini	2,22	45°	10,18	2,67	5,41	2,10	M1,6	20 Ncm
	19.273	TIN screw hex Standard	2,33	45°	8,35	2,60	2,95	2,80	M2,0	30 Ncm
	19.672	TIN screw hex Mini	2,22	45°	10,18	2,67	5,41	2,10	M1,6	20 Ncm
	19.673	TIN screw hex Standard	2,33	45°	8,35	2,60	2,95	2,80	M2,0	30 Ncm




INTERNAL BT

compatible with
BTI®

SOCKET
HEX 1,22mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.065	Screw hex MU	2,00	0°	3,62	1,50	0,32	1,80	M1,4	10 Ncm
	19.079	Screw hex	2,55	0°	7,05	1,80	2,05	3,20	M1,8	35 Ncm
	19.665	TIN screw	2,00	0°	3,62	1,50	0,32	1,80	M1,4	10 Ncm
	19.679	TIN screw MU	2,55	30°	7,05	1,80	2,05	3,20	M1,8	35 Ncm



A. HEAD DIAMETER (mm)
 B. HEAD BASE ANGLE
 C. OVERALL LENGTH (mm)
 D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
 F. THREAD LENGTH (mm)
 G. THREAD KIND

**INTERNAL SEM**compatible with
Bego Semados®SOCKET
HEX 1,27mm**DESS**

> DENTAL SMART SOLUTIONS




REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.091	Screw hex 3.25/5.5	2,41	55°	7,65	2,47	1,38	3,80	M1,8	30 Ncm
19.691	TIN screw hex 3.25/5.5	2,41	55°	7,65	2,47	1,38	3,80	M1,8	30 Ncm

**MULTI-UNIT
CONNECTION**SOCKET
HEX 1,20mmSOCKET
HEX 1,27mmSOCKET
UniGrip®

REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.006	Screw UniGrip® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
19.031	Screw UniGrip® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 Ncm
19.055	Screw hex 1,20	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
19.065	Screw hex 1,22	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
19.085	Screw hex Multi-Unit 4.8	2,15	0°	3,50	1,30	0,00	2,20	M1,6	25 Ncm
19.110	Screw hex 1,27 RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
19.206	DLC screw UniGrip® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
19.231	DLC screw UniGrip® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 Ncm
19.606	TIN screw hex 1,20 on 3i® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
19.631	TIN screw UniGrip® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 Ncm
19.685	TIN screw hex Multi-Unit 4.8	2,15	0°	3,50	1,30	0,00	2,20	M1,6	25 Ncm

Screw C-Base®SOCKET
HEX 1,00mmSOCKET
HEX 1,27mm

REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.052	Screw hex 1,00 C-Base®	2,15	45°	8,87	1,91	4,16	2,80	M1,6	20 Ncm
19.702	Screw hex 1,27 C-Base®	1,90	60°	9,37	1,50	5,97	3,40	M1,8	15 Ncm



A. HEAD DIAMETER (mm)

B. HEAD BASE ANGLE

C. OVERALL LENGTH (mm)

D. HEAD LENGTH (mm)

E. DISTANCE FROM
HEAD TO THREAD (mm)

F. THREAD LENGTH (mm)


G. THREAD KIND



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.306	Screw MUA AURUMBase®	2,20	0°	3,60	1,50	0,10	2,00	M1,4	15 Ncm
	19.331	Screw UNI AURUMBase®	2,30	0°	4,10	1,40	0,70	2,00	M1,8	15 Ncm
	19.400	Screw AURUMBase®	2,30	0°	7,40	1,40	3,30	2,70	UNF1,72	30 Ncm
	19.401	Screw AURUMBase®	2,10	0°	6,10	1,40	2,50	2,20	M1,6	35 Ncm
	19.402	Screw AURUMBase®	2,50	0°	5,90	1,40	1,69	2,81	M2,0	35 Ncm
	19.403	Screw AURUMBase®	2,30	0°	5,90	1,40	2,50	2,00	M1,8	35 Ncm
	19.404	Screw AURUMBase®	2,30	0°	7,20	1,40	3,55	2,25	M1,8	35 Ncm
	19.405	Screw AURUMBase®	2,45	0°	8,05	1,40	4,60	2,05	M2,0	35 Ncm
	19.417	Screw AURUMBase®	2,30	0°	4,90	1,40	2,90	3,85	UNF-1-72	30 Ncm
	19.425	Screw AURUMBase®	2,45	0°	8,90	1,40	5,45	2,00	M2,0	25 Ncm
	19.438	Screw AURUMBase®	2,10	0°	7,90	1,40	4,50	2,05	M1,6	● Ncm
	19.441	Screw AURUMBase®	2,10	0°	7,30	1,40	3,70	2,00	M1,6	35 Ncm
	19.442	Screw AURUMBase®	2,45	0°	7,20	1,40	3,75	2,05	M2,0	35 Ncm
	19.443	Screw AURUMBase®	2,10	0°	8,15	1,40	4,75	2,00	M1,6	35 Ncm
	19.444	Screw AURUMBase®	2,10	0°	9,20	1,40	5,60	2,20	M1,6	● Ncm
	19.446	Screw AURUMBase®	2,45	0°	6,70	1,40	3,25	2,05	M2,0	35 Ncm
	19.460	Screw AURUMBase®	2,30	0°	7,42	1,40	4,02	2,00	M1,8	25 Ncm
	19.461	Screw AURUMBase®	2,45	0°	7,75	1,40	4,30	2,05	M2,0	25 Ncm
	19.468	Screw AURUMBase®	2,10	0°	6,90	1,40	3,05	2,45	M1,4	15 Ncm
	19.474	Screw AURUMBase®	2,10	0°	8,70	1,40	5,30	2,00	M1,6	20 Ncm
	19.476	Screw AURUMBase®	1,90	0°	8,50	1,40	5,05	2,05	M2,0	20 Ncm

● Internal Hex "Click"= 25 Ncm
Internal Hex FD= 24 Ncm

● Internal CAM= 20 Ncm
Internal OSS (Mini)= 24 Ncm
Internal ICX= 30 Ncm



A. HEAD DIAMETER (mm)

B. HEAD BASE ANGLE

C. OVERALL LENGTH (mm)

D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)

F. THREAD LENGTH (mm)

G. THREAD KIND



www.dessdental.com