

3M[™] Scotchbond[™] Universal Plus Adhesive

The original, now better than ever!



The benchmark adhesive you trust.

3M[™] Scotchbond[™] Universal Plus Adhesive takes universal adhesives to the next level. It keeps all the benefits of the original 3M[™] Scotchbond[™] Universal Adhesive while adding even greater control and predictability – because a good adhesive should provide more than just a strong bond.

S steps forward from the original

Added radiopacity

Bonds and seals caries-affected dentine

Advanced bonding to dental and restorative substrates, including glass ceramics

Full dual- and self-cure compatibility without additional dual-cure activator bottle

Free of BPA derivatives like BisGMA



3M Scotchb

niversal Plus dhesive



3 Neuss-C

Plus, the original Scotchbond Universal Adhesive's popular features:

Universally suitable

- 🕑 Direct
- S Indirect
- All etching techniques

. now setting a new standard


Minimally Invasive Dentistry

Added radiopacity

Radiolucency under an existing restoration could lead to misdiagnosis and invasive overtreatment. 3M[™] Scotchbond[™] Universal Plus Adhesive is the first radiopaque universal adhesive with dentine-like radiopacity, helping to overcome this issue.

These doubtful X-rays happen every day; most often critical step in bitewing radiograph, class II in proximal box."

- Prof. Dr. Angelo Putignano, University of Marche-Ancona, Italy



Intentional adhesive pooling under Class I restorations - X-rays of extracted human molars, filling material 3M™ Filtek™ One Bulk Fill Restorative. K. Dede et al., J Dent Res 99 (Spec Iss A): 192 IADR, 2020 - 3M internal data, H. Loll et al., J Dent Res 99 (Spec Iss A): No. 757, 2020.


Radiolucent adhesive pooling underneath a restoration doubles the risk of restoration replacement."

- Dr. John O. Burgess, University of Alabama, USA


Bonds and seals caries-affected dentine to help preserve tooth structure

Scotchbond Universal Plus Adhesive gives you confidence when performing minimally invasive preparations because it has the same high bond strength to cariesaffected dentine as to sound dentine.

In addition to the excellent quality of the restoration margins, it seals caries-affected dentine by forming a well-defined, void-free hybrid layer - because "the seal is the deal" to prevent secondary decay.



Shear bond strength on caries-affected and sound dentine - C. Thalacker, et al., Interaction of an Experimental Universal Adhesive with Caries-Affected Dentin, J Dent Res 99 (Spec Iss A): 191 IADR, 2020.



Caries zones - Image courtesy of Prof. Leandro Hilgert and Prof. Soraya Leal, University of Brasilia (UnB), Brazil.

Bonds all restorative materials


without the need for separate primers

Virtually no post-operative sensitivity

of dentists reported that 3M[™] Scotchbond[™] **9%** Universal Plus Adhesive shows virtually no post-operative sensitivity*

Advanced bonding to glass ceramic

Scotchbond Universal Plus Adhesive reaches the gold standard level of traditional, separate silanes.



Report conclusion:

3M[™] Scotchbond[™] Universal Plus Adhesive performed better than 3M[™] Scotchbond[™] Universal Adhesive for glass ceramic (IPS e.max[®] CAD, Ivoclar Vivadent) bonding. After the glass ceramic was etched with HF acid, a favorable 'immediate' and 'aged' share bond strength were recorded.

Scotch

Shear bond strength on glass ceramic - Prof. Dr. Bart van Meerbeek et al., KU Leuven, Belgium, Report to 3M, 2020

Fully aligned system with 3M™ RelyX™ Universal Resin Cement

Excellent bond strength for virtually all dual-cure resin cement indications.[^] Simplified workflows – 3M Scotchbond Universal Adhesive is cured when in contact with RelyX Universal Resin Cement.

Ordering Information


Item no. Product information

41293	Intro Kit Vial 1 Vial 5ml, 50 Disposable Brushes, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41294	Refill Vial – 1 x 5ml
41295	Refill Vial – 3 x 5ml
41263	3M[™] Scotchbond[™] Universal Etchant Etching Gel 2 Syringes x 3ml, 50 Dispensing Tips

41297	Intro Kit Unit Dose 50 Unit Doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41298	Unit Dose 100 100 Unit Doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41299	Unit Dose 200 200 Unit Doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips



Vial with smooth and even surface - easy to wipe off excess



3M.com.au/ScotchbondPlus



Scan to find out more

3M Australia Pty Limited Building A, 1 Rivett Rd North Ryde, NSW 2113 1300 363 454 www.3M.com.au **3M New Zealand Limited** 94 Apollo Drive Rosedale, Auckland 0632 0800 808 182 www.3M.co.nz

^Resin Cement Bond Strength to Multiple Substrates. Dental Advisor, M. Cowen, J.M. Powers, 2020. 3M, Filtek, RelyX and Scotchbond are trademarks of 3M Company. © 3M 2021. All rights reserved. Please recycle. AV011477542.