Science. Applied to Life.™

3M[™] Scotchbond[™] Universal Plus Adhesive

More than just a strong bond.



3M™
Scotchbond™
Universal Plus
Adhesive takes
universal adhesives
to the next level.
It keeps all the benefits of the
original 3M™ Scotchbond™
Universal Adhesive, the first truly
universal adhesive, while adding
even greater control and
predictability – because
a good adhesive should
provide more than just
a strong bond.

From excellent...

Built on the clinically proven and world's most researched universal adhesive.

- 0
- All direct and indirect indications light-curing, self-curing or dual-curing methacrylate-based composites / cements
- All etching techniques
 total-etch, self-etch, and selective
 enamel etch
- Bonds to all restorative materials without the need for separate primers
- Available in vial and unit dose for efficient hygiene management
- Virtually no post-operative sensitivity in field evaluation:

99% of dentists reported that 3M™ Scotchbond™ Universal Plus Adhesive shows virtually no postoperative sensitivity.*

... to exceptional

Formulation carefully optimised to achieve unique features.



First radiopaque universal adhesive reducing the risk of X-ray misdiagnosis and invasive overtreatment



Bonds and seals caries-affected dentine maximising preservation of natural tooth structure, supporting minimally invasive preparation guidelines



Enhanced Silane Functionality that offers Gold Standard adhesion to all dental substrates and in particular, glass ceramics



Full dual- and self-cure compatibility no need for an additional dual-cure activator bottle



BPA free formulation free of BPA derivatives like BisGMA






Fully aligned system with 3M™ RelyX™ Universal Resin Cement

excellent bond strength for virtually all resin cement indications. 3M Scotchbond Universal Plus Adhesive is cured when in contact with RelyX Universal Resin Cement



Added radiopacity

On average, dentists see at least one difficult to interpret X-ray per day due to radiolucency under an existing restoration. This could lead to misdiagnosis and invasive overtreatment. 3M™ Scotchbond™ Universal Plus Adhesive is the first radiopaque universal adhesive with dentine-like radiopacity, helping to overcome this issue.



Intentional adhesive pooling under Class I restorations - X-rays of extracted human molars, filling material 3M™ Filtek™ One Bulk Fill Restorative. 3M internal data, H. Loll et al., J Dent Res 99 (Spec Iss A): No. 757, 2020

Radiolucent adhesive pooling underneath a restoration doubles the risk of restoration replacement."

Dr. John O. Burgess, University of Alabama, USA

Ordering information

Item no. Product information

41293 Intro Kit Vial 1 Vial 5ml, 50 Disposable Brushes, 1 Etchant Syringe 3ml, 25 Dispensing Tips

41294 Refill Vial - 1x 5ml

Refill Vial 3 x 5ml 41295

41263 3M™ Scotchbond™ Universal Etchant Etching Gel

2 Syringes x 3ml, 50 Dispensing Tips

3M™ Scotchbond™ Universal Plus Adhesive - Unit Dose

41297 **Intro Kit Unit Dose** 50 Unit Doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips

41298 Unit Dose 100

100 Unit Doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips

41299 Unit Dose 200

200 Unit Doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips



with smooth and even surface easy to wipe off excess



for efficient hygiene management reduces the risk of cross contamination

3M.com.au/ScotchbondPlus







Scan to find