Science. Applied to Life.™

3M™ RelyX™ Universal Resin Cement & 3M™ Scotchbond™ Universal Plus Adhesive

The two that can do it all.

3M™ RelyX™ Universal Resin Cement works both as a standalone, self-adhesive cement - and as an adhesive cement together with 3M™ Scotchbond™ Universal Plus Adhesive. The adhesive serves as a universal primer for all materials and covers all direct and indirect bonding indications.

Eliminate the hassle and confusion of multiple resin cements, primers and adhesives and simplify your direct and indirect restorative workflows. Less products to keep on stock and replenish. Saves space, time and money.

- First radiopaque all-inone universal adhesive
- Further boosts bond strength of RelyX **Universal Resin** Cement
- Fully aligned system: adhesive cured by cement, no light-cure needed
- **Universal primer** for all restorative materials



Game-changing syringe. Proprietary 3M design.

- Micro Mixing Tip providing 80% less cement waste*
- 15 20 applications from 3.4g**
- Hygienic self-sealing syringe - for storage without the used mixing tip.



- More ergonomic shape with 50% less plastic waste*
- Thin, long elongation tip for easy root canal









Excellent self-adhesive bond strength to dentine in self- and light-cure

Shear bond strength (MPa) to dentin after artificial aging





Ordering information

3M™ RelyX™ Universal Resin Cement

Item no. Product information

56968 4 x 3M™ RelyX™ Universal Resin Cement syringes shades TR, A1, A3O, WO (3.4g each), 60 x 3M™ RelyX™ Universal Micro Mixing Tips, 20 x Elongation Tips, 1x vial 3M™ Scotchbond™ Universal Plus Adhesive (5 ml), 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 x orange microbrushes, step-by-step guides

Trial Kits

Each trial kit contains: 1 x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15 x 3M™ RelyX™ Universal Micro Mixing Tips, 5 x Elongation Tips, 1 x vial 3M™ Scotchbond™ Universal Plus Adhesive (1.5 ml), 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 x blue etchant tips, 50 x orange microbrushes, step-by-step guides

56969 Shade TR

56970

Shade A1

Each refill contains: 1 x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15 x 3M™ RelyX™ Universal Micro Mixing Tips, 5 x Elongation Tips

56971 Shade TR 56973

Shade A3O

56972 Shade A1 56974 Shade WO

Value Packs

Each value pack contains: 3 x 3M™ RelyX™ Universal Resin Cement syringes (3.4g each), 45 x 3M™ RelyX™ Universal Micro Mixing Tips, 15 x Elongation Tips

56977 Shade TR 56978

Shade A1

Mixing Tips

30 x 3M™ RelyX™ Universal Micro Mixing Tips 56975

15 x 3M™ RelyX™ Universal Micro Mixing Tips with Elongation Tips 56976

3M™ Scotchbond™ Universal Plus Adhesive

41293 Intro Kit: 1 x 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml), 50 x disposable brushes, 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 x dispensing tips

Refill: 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml)

1x5ml vial 41294

41295

3 x 5 ml vial

Unit Dose

3M™ Scotchbond™ Universal Plus Adhesive Unit Dose L-Pops, 1 x 3M[™] Scotchbond[™] Universal Etchant (3 ml), 25 x dispensing tips

41297 Intro Kit: 50 x L-Pop 41299

200 x L-Pop

41298 100 x L-Pop

3M™ Scotchbond™ Universal Etchant Etching Gel

2 x 3M™ Scotchbond™ Universal Etchant (3 ml each), 50 x dispensing tips



3M™ RelyX™ Universal Resin Cement Intro Kit (56968)



3M™ RelyX™ Universal Resin Cement Trial Kit TR (56969)



3M™ RelvX™ Universal Resin Cement Refill TR (56971)



3M™ Scotchbond™ Universal Plus Adhesive Refill (41294)



3M™ Scotchbond™ Universal Plus Adhesive **Unit Dose:**

For efficient hygiene management. Reduces the risk of cross-contamination.

3M™ Scotchbond™ Universal Plus Adhesive Unit Dose 100 L-Pop (41298)