

SOFTLY

Universal mixer for predosed capsules

Softly is a universal high-speed mixer for predosed dental materials in capsules (amalgams, glass ionomers).

It is controlled by a microprocessor which means precise repeatability of the selected time to guarantee you a controlled, consistent mix.

The esthetically-pleasing ergonomic design of the **Softly** allows for easy use and selection of the timer display.

The width of the oscillating range (25mm), the mixing frequency (4 000 vibrations/min ± 50 /minute) and your selection of the correct time required, provide perfect homogeneity and an ideal consistency for the type of filling material that you choose.

Softly is very stable, extremely quiet, and very easy to clean and maintain.

To guarantee totally safe operation, **Softly** is deactivated when you open the front cover to access the capsule.

Softly is made entirely of injection-molded ABS material. It complies with all electrical and mechanical safety standards.



Precise repeatability of the selected time



Deactivated when the front cover is open



Ideal consistency whatever the type of dental material you choose



Features



High mixing frequency
Noiseless operation
Excellent stability
All-electric safety mechanism
Ergonomic design
Maintenance-free
Easy to clean

TECHNICAL DATA

Supply voltage (upon request):	230V~ - 50Hz
Absorbed power:	
Fuse:	
Mixing frequency:	,
Amplitude:	
Operating conditions:	intermittent duty
Programmable mixing time:	•
Overall dimensions (W x D x H):	
Weight.	3 1kg



www.acteongroup.com

