

# SDI | RADIOMETER X

LED LIGHT CURING METER




## OUTSTANDING ACCURACY. NO CALIBRATION. NO DELAY.

Compatible with all LED curing lights, the new Radiometer X accurately checks LED intensity to confirm your curing performance.

- **No calibration:** Simply point and measure. No need to enter dimensions
- **No dependence:** No pairing, No apps, No display delay
- 5 year warranty

## STUDY: PROVEN ACCURACY

Results from a study of 7 radiometers and 23 curing lights has confirmed the outstanding accuracy of the Radiometer X.



Source: Dr Carlos Jose Soares [PhD]

~ Radii Plus; Radii Cal; Radii Xpert/Monowave tip; Radii Xpert/Multiwave tip.

\*\* Bluephase Meter II includes a LED diameter checker. Accuracy readings will vary if the user enters an incorrect diameter or includes the LED casing as part of the diameter.

COMPOSITES | GLASS IONOMERS | TOOTH WHITENING | SDF RIVA STAR | CEMENTS | ADHESIVES | ETCHANTS | SEALANTS | AMALGAMS | ACCESSORIES | EQUIPMENT