



# Fine Gauge, Short Needle

AVAILABLE  
NOW

For more comfortable injections

- 31G ultra-thin wall needle
- Ultra-sharp double bevel for easy skin penetration
- Easy to read markings on the clear barrel promotes accurate dosing
- Low dead space design minimises waste
- Latex and DEHP free
- Individually wrapped
- Sterile



## 0.3mL Fixed Needle Syringe

Product Code	Gauge	Volume		Length mm	Graduations		Inner Box Qty Each	Carton Qty Each
		mL	Units		mL	Units		
U31PK	31G	0.3	30	8	0.01	1	100	3600

Distributed by Terumo Australia Pty Ltd  
Macquarie Park NSW 2113 Australia T: 1800 837 866  
Mangere Auckland 2022 New Zealand T: 0800 66 77 57

All products manufactured for Exchange Supplies, UK  
GHP220603D1439V2 © 2022 Terumo Australia Pty Ltd

