



What doctors are saying:

"In a single day I was able to do 11 mock-up veneers that were not part of the original schedule, resulting in several new cases being accepted. I wouldn't want to work without them now."

-Dr. Chad Wagener, DDS

"Terrific tool to quickly and easily create beautiful anterior restorations."

-Dr. Gary M. Radz, DDS

"The simplicity of Uveneer is absolutely remarkable. Why didn't I think of this?"

-Dr. John C. Comisi, DDS, MAGD

"Uveneer makes the practice of the anterior esthetic dentistry easier, faster, and better."

-Dr. George Freedman, DDS





UVKEV1 Uveneer Extra Kit

- 6 x Extra Large upper anterior templates 6 x Large upper anterior templates
- 6 x Medium upper anterior templates
- 6 x Square upper anterior templates

UVKEXLSQV1 Uveneer Extra XL & SQ Kit

6 x Extra Large upper anterior templates

6 x Square upper anterior templates

UVKELMV1 Uveneer Extra L & M Kit

6 x Large upper anterior templates 6 x Medium upper anterior templates

Follow Us!

Join the Uveneer template Facebook group to share cases, learn techniques, and ask Dr. Sigal Jacobson questions! f facebook.com/groups/uveneer

Follow the Uveneer template brand for news, updates, and videos.

f facebook.com/uveneerdentist









DIRECT COMPOSITE TEMPLATE SYSTEMS



Simple and Versatile

The Complete Uveneer Template Family

The Uveneer and Uveneer Extra™ template systems are unique, minimally invasive, and create beautiful direct composite veneers with predictable shape and symmetry. They mimic natural dentition and are designed to create high-quality anterior restorations in one visit. They can also be used for cosmetic mock-ups and shade selection, as well as temporaries during porcelain veneer creation.



Uveneer Templates Promise a Beautiful Symmetry



A Way to Provide Unique and Natural Smiles



The original Uveneer template kit contains upper and lower anterior templates in Medium and Large which:

- Create predictable, reproducible composite restorations
- · Prevent the oxygen inhibition layer during curing, resulting in a hard and glossy surface
- Allow light to pass through the template to the composite for effective curing
- · Work with any preferred composite
- Release easily from cured composite resin
- Facilitate application on individual or multiple teeth
- Are autoclavable* and reusable, making them a cost-effective choice

*Black base and transparent lids are NOT autoclavable. Only templates are autoclavable

Uveneer Extra templates provide the option of four new upper anterior designs: Extra Large, Large, Medium, and Square. To complement the original Uveneer templates, they also:

- Create natural-looking contours on composite restorations
- Require minimal adjusting or polishing, saving time
- Mimic real tooth anatomy for authentic, natural results
- Are molded based on scanned teeth from Dr. Jan Hajto's work
- Add to the unique potential of the patient's smile
- · Possess all of the trusted qualities of the original Uveneer templates, but with a new design





















L - Large Upper

L - Large Lower

M - Medium Upper M - Medium Lower

XL - Extra Large

L - Large

M - Medium

SQ - Square