

SANIFLEX PTY LTD

## PRODUCT SHEET

### NITRILE EXAMINATION GLOVES

100 PACK

PRODUCT IMAGE	DESCRIPTION
	<p><b>NITRILE EXAMINATION GLOVES - POWDER FREE</b>  ITEM NO.: NIT100B  UNITS PER BOX: 100 GLOVES PER BOX  BOX DIMENSIONS: 20 x 11 x 6cm  BOXES PER CARTON: 10  UNITS PER CARTON: 1000  CARTON DIMENSIONS: 31.5 x 23.5 x 21cm  SHELF LIFE: 5 Years  STORAGE INSTRUCTIONS: Keep out of direct sunlight, store in a cool and dry place. Keep away from sources of ozone and ignition.</p>
	<p><b>PRODUCT INFORMATION</b>  MATERIAL: 100% Nitrile  SIZES: XS, S, M, L, XL  COLOUR: Blue  GLOVE LENGTH: ≥240mm  FINGER THICKNESS(mm): ≥0.05  PALM THICKNESS (mm): ≥0.05  FORCE AT BREAK: ≥6N before and after aging  TENSILE STRENGTH: ≥14MPa before and after aging  ELONGATION: 500% minimum before aging, et 400% after aging  Dimensions(±10mm) XS/70;S/80;M/95;L/110;XL/120  SURFACE: Micro-textured on fingers</p>
	<p><b>INTENDED USE:</b>  Saniflex nitrile examination gloves are intended to be worn on the hands of healthcare personnel to prevent contamination between patients and examiners.</p>
	<p><b>FEATURES:</b>  Manufactured under ISO 9001:2015 and ISO 13485:2016 Quality Management System  Textured fingertips for enhanced grip  Custom design enhanced comfort and fit  An alternative solution for individuals who are allergic to natural rubber latex</p>
	<p><b>STANDARDS &amp; CERTIFICATIONS:</b>  COMPLIANT WITH EN455-1:2020, EN455-2:2015, EN455-3:2015, EN455-4:2019  TGA, HACCP</p>