

SANIFLEX PTY LTD

PRODUCT SHEET

Harley N95 Surgical Respirator with Fluid Resistance

PRODUCT IMAGE	DESCRIPTION
 <p>The image shows a box of Harley N95 Surgical Respirators and a single mask. The box is white with a red top section and features the text '20 Masks', 'N95 Particulate Respirator', and 'MODEL NO L-188'. The mask is white with yellow ear hooks and a yellow adjustable nose clip. The mask has 'Harley N95 L-188 Surgical Respirator' printed on it.</p>	<p>HARLEY N95 SURGICAL RESPIRATOR WITH FLUID RESISTANCE, 20 masks per box</p> <p>ITEM NO.: N95-H20</p> <p>UNITS PER CARTON: 20</p> <p>Material: PP Non-woven fabric: High efficiency low-resistance melt blown, electrostatic cotton.</p> <p>Headbands: Latex free</p> <p>Specifications:</p> <ul style="list-style-type: none"> At least 99% Particle Filtration Efficiency Adjustable nose clip Five dimensional Thickened material Comfortable ear hook <p>ARTG NO.: 341498</p>