



MANAGING PAIN FOR YOUR PRACTICE



SEPTANEST 1:100,000

4% articaine with 1:100,000 epinephrine injection solution

Indications

Local or loco-regional dental anaesthesia in patients of at least 4 years old

Features and Benefits

- Septanest 40 mg/ml plus adrenaline 0.01 mg/ml (1/100 000) is particularly appropriate for procedures of longer duration and when there is a risk of significant bleeding into the operative field
- Fast onset to avoid delay before effective analgesia is achieved
- No-methylparaben formulation to avoid allergic reactions on sensitive patients
- 100% latex-free
- Terminal sterilization of product/cartridges, the Gold Standard method for Health Authorities around the world
- Mylar wrapped cartridges to minimize risk of injury in the event of a cartridge breakage with blister packaging to avoid contamination

Presentation

- 2 x Box of 50 2.2ml cartridges



LIGNOSPAN SPECIAL

2% lidocaine with 1:80,000 epinephrine injection solution

Indications

Local or loco-regional dental anaesthesia

Features and Benefits

- No-methylparaben formulation to avoid allergic reactions on sensitive patients
- 100% latex-free
- Terminal sterilization of product/cartridges, the Gold Standard method for Health Authorities around the world.
- Mylar wrapped cartridges to minimize risk of injury in the event of a cartridge breakage with blister packaging to avoid contamination

Presentation

- 2 x Box of 50 2.2ml cartridges



SEPTANEST 1:200,000

4% articaine with 1:200,000 epinephrine injection solution

Indications

Local or loco-regional dental anaesthesia in patients of at least 4 years old

Features and Benefits

- Septanest containing 1/200 000 epinephrine should be preferred as it reduces the risk of excessive uptake of Septanest into the systemic circulation
- No-methylparaben formulation to avoid allergic reactions on sensitive patients
- 100% latex-free
- Terminal sterilization of product/cartridges, the Gold Standard method for Health Authorities around the world
- Mylar wrapped cartridges to minimize risk of injury in the event of a cartridge breakage with blister packaging to avoid contamination

Presentation

- 2 x Box of 50 2.2ml cartridges



SCANDONEST 3% PLAIN

3% mepivacaine injection solution

Indications

Local or loco-regional dental anaesthesia

Features and Benefits

- Formulation without vaso-constrictor for specific medical situations
- No-methylparaben formulation to avoid allergic reactions on sensitive patients
- Terminal sterilization of product/cartridges, the Gold Standard method for Health Authorities around the world
- 100% latex-free
- Mylar wrapped cartridges to minimize risk of injury in the event of a cartridge breakage with blister packaging to avoid contamination

Presentation

- 2 x Box of 50 2.2ml cartridges



SEPTOJECT XL NEEDLES

Sterile single use dental needle with enlarged bore

Indications


Routine administration of dental local anaesthetics

Features and Benefits

- With enlarged bore (43% wider than the one of a standard needle) to reduce level of pressure during injection, thus leading to less pain for the patient
- Less trauma to the tissue achieved through a triple-bevel needle, a polished & siliconised canula and an excellent flexibility of the needle
- High grade surgical stainless steel tubing to reduce risk of breakage
- With bevel indicator to assist in orientating the bevel during injection

Presentation

- Box of 100 Sterile Needles



Length	Gauge	Colour Code
X Short - 10 mm	30 G	Purple
Short - 25 mm	30 G	Blue
Short - 25 mm	27 G	Red
Long - 35 mm	27 G	Yellow

EVOLUTION NEEDLES

The innovative needle with a patented scalpel-designed bevel

Indications

Routine administration of local dental anesthetics by infiltration or intraligamentary injections

Features and Benefits

- Engineered for a smoother penetration with less tissue displacement, less discomfort for patients even with multiple injections
- Less force required combined with significantly reduced deflection, better control, and accuracy

Presentation

- Box of 100 Sterile Needles



Length	Gauge	Colour Code
X Short - 9 mm	30 G	Purple
Short - 25 mm	30 G	Blue
Short - 25 mm	27 G	Red

SEPTOJECT NEEDLES

Indications

Routine administration of dental local anaesthetics

Features and Benefits

- Less trauma to the tissue achieved through a triple-bevel needle, a polished & siliconised canula and an excellent flexibility of the needle
- High grade surgical stainless steel tubing to reduce risk of breakage

Presentation

- Box of 100 Sterile Needles



Length	Gauge	Colour Code
X Short - 16 mm	30 G	Purple
Short - 25 mm	30 G	Blue
Short - 25 mm	27 G	Red
Long - 35 mm	27 G	Yellow

SEPTOJECT SELF ASPIRATING SYRINGE

Indications

Routine administration of dental local anaesthetics

Features and Benefits

- Made from stainless steel, these are suitable for any sterilisation
- Designed especially for 2.2ml cartridges with a hole in the plunger
- Made in one piece, no parts to lose
- Light weight to reduce hand fatigue during injections



ULTRA SAFETY PLUS TWIST

Indications

Needle stick injury prevention device

Features and Benefits

- Protection from needle stick injuries for you and your staff
- Easy and secure lock makes your device intuitive and safe
- Both passive and active aspiration for the security of the patient
- Transparent barrel to allow visualization of aspiration
- Bevel indicator facilitating the appropriate orientation of the bevel
- Two different positions: Holding Position (reversible) and locking position (irreversible)
- Available with either sterile WHITE single-use (SP-10698Q) or sterilizable BLUE handles (SP-10697P)
- Complies with latest regulations

Presentation

- Box of 50 Sterile Single use (white) Handles
- Single Sterilizable (Blue) autoclavable handle with 5 O-rings
- Box of 100 Sterile Needles + 1 Sterilizable (Blue) Handle

Length	Gauge	Colour Code
X Short - 10 mm	30 G	Purple
Short - 25 mm	30 G	Blue
Short - 25 mm	27 G	Orange
Long - 35 mm	27 G	Yellow



Indications

- Topical analgesia of oral mucosa before anaesthetic injection
- Temporary pain relief during short and simple dental procedures
- Gag reflex elimination when performing radiographs or impressions

XOGEL ENFANT

5% lidocaine analgesic gel

- Gel consistency prevents wash away to improve activity
- Improves patient's comfort during many procedures
- Cherry flavour

Presentation

- 50g dispenser bottle



XYLONOR GEL

5% lidocaine analgesic gel

- Gel consistency prevents wash away to improve activity
- Fast acting, gives anaesthesia within 2-5 mins
- Mint flavoured gel

Presentation

- 15g tube



XYLONOR PELLETS

- Rapid acting, producing anaesthesia in 10-15 seconds and deep anaesthesia in 30-60 seconds
- Apply to area with tweezer using painting action
- Gives rapid anaesthesia within 2-5minutes

Presentation

- Jar of 200 Lignocaine base cotton pellets





MANAGING PAIN FOR YOUR PRACTICE

