

LIMITED TIME ONLY - WHILE STOCKS LAST

# Simplicity

by **Kerr™**

3 solutions to simplify your complete restorative workflow

## SimpliShade™

This universal composite covers all 16 VITA® classical shades with just 3 shades, and is powered by Kerr's Adaptive Response Technology. No blocker required!




## OptiBond™ Universal

One-bottle universal adhesive, providing excellent adhesion to a variety of substrates for direct and indirect applications - even indirect metal-based restorations when used with Maxcem Elite™.




## Maxcem Elite™

This cement now features universal flexibility as a self-adhesive cement, or as an adhesive-resin cement when paired with OptiBond™ Universal on a preparation requiring more retention.



**NOW INCLUDED**

**FREE**  
**OptiDisc® Mini Kit** (OptiDisc Mini Kit REF/Part No. 4188)



### Simplicity Kit contains:

- 5x SimpliShade™ Syringes
- 1x OptiBond™ Universal
- 1x Maxcem Elite™

(Simplicity Kit REF/Part No. 377777)

**Save \$300 on RRP** **Kerr™**

# Ordering Information

## SimpliShade™ Universal Composite




<b>37001</b>	Syringe, Light	<b>37006</b>	Unidose 20 Pack, Light
<b>37002</b>	Syringe, Medium	<b>37007</b>	Unidose 20 Pack, Medium
<b>37003</b>	Syringe, Dark	<b>37008</b>	Unidose 10 Pack, Dark
<b>37004</b>	Syringe, Bleach White	<b>37009</b>	Unidose 10 Pack, Bleach White
<b>37005</b>	Syringe, Universal Opaque	<b>37010</b>	Unidose 10 Pack, Universal Opaque

## OptiBond™ Universal Bonding Agent



<b>36517</b>	Universal Bottle Kit
<b>36518</b>	Universal Unidose (100 Pack)
<b>36519</b>	Universal Bottle Refill
<b>26909</b>	Gel Etchant Dispensing Tips (pkg of 30)
<b>31297</b>	Gel Etchant Refill Pack

## Maxcem Elite™ Universal Resin Cement



<b>34417</b>	Maxcem Elite Intro Kit
<b>33872 / 34418</b>	Clear (2 / 4 dual syringes)
<b>33873</b>	White (2 dual syringes)
<b>33874</b>	White Opaque (2 dual syringes)
<b>33875</b>	Yellow (2 dual syringes)
<b>32862</b>	50 4:1 automix tips (regular)
<b>33879</b>	8 4:1 automix tips (wide)
<b>33880</b>	8 intraoral tips
<b>33881</b>	8 root-canal tips

## OptiDisc™ Mini Kit



<b>4188</b>	OptiDisc Mini Kit ø 12.6 mm: 30 x Extra-coarse, 30 x Coarse, 30 x Fine, 30 x Extra-fine, 1 standard mandrel, 1 short mandrel, 1 OptiShine
-------------	---

**Scan the QR code to order your Simplicity Kit  
with bonus OptiDisc® Mini Kit.  
Your Kerr Representative will contact you.**

