

# TRIOS 3 Wired

Proven. The trusted choice.

- ✓ High accuracy and precision
- ✓ Invisalign connected
- ✓ Realistic scan colors
- ✓ Automatic scan clean-up
- ✓ Print-ready model design
- ✓ Unrivalled open dental-ecosystem
- ✓ Patient apps to boost case acceptance
- ✓ Opportunity for in-house design and production

## TECHNICAL SPECIFICATIONS

### GENERAL

<b>What's in the box</b>	TRIOS 3 scanner wand, TRIOS 3 scanner pod, 5x TRIOS 3 scanner tips, TRIOS 3 protection tip, TRIOS 3 calibration tip, TRIOS 3 color calibration tip, TRIOS 3 color calibration kit, 3Shape USB-C custom cable, 3Shape license dongle, USB dongle with software, documentation
<b>Packaging Dimensions (LxWxH cm)</b>	44 x 12,5 x 31
<b>Invisalign Connection</b>	Yes (except NA, China and Japan)

### HARDWARE

<b>Scanner Weight</b>	340 g
<b>Scanner Dimensions (LxWxH cm):</b>	27,3 x 4 x 4,8
<b>Tip Dimensions (LxWxH cm):</b>	8,6 x 2 x 2,05
<b>Scanner Cord Length (cm)</b>	150
<b>Pod Dimensions (DxH cm)</b>	12 x 7
<b>Pod Weight</b>	510 g
<b>USB-C Cable Length (cm)</b>	100



Let's change dentistry together