

# TRIOS 6 Wireless

Revolutionary. TRIOS Redefined.

- ✓ Clinically proven accuracy and precision
- ✓ Best-ever TRIOS scan quality, with up to 110% higher scan resolution
- ✓ Hyperspectral technology capturing white, fluorescent, and near-infrared data simultaneously in one scan
- ✓ TRIOS 6 scans are compatible with TRIOS Dx Plus for AI assistive detection of most common oral health conditions<sup>1</sup>
- ✓ Empowered patient communication with DentalHealth app
- ✓ Full mobility with wireless and calibration-free scanning
- ✓ Optimal-power battery charging
- ✓ Energy-saving sleep mode
- ✓ ScanAssist intelligent-alignment technology
- ✓ Sensory feedback guiding scanning
- ✓ Next-level ergonomics and hygienic design, for minimized patient cross contamination
- ✓ Unrivalled open dental-ecosystem and patient apps to boost case acceptance

## Technical specifications

### GENERAL

|                   |  |
|-------------------|--|
| What's in the box | TRIOS 6 scanner, 1 autoclavable tip, a box of 25 TRIOS Ready Tips, TRIOS protection tip, a box of TRIOS Sleeves, 3 batteries, charger including power cord, scanner pod, wireless adapter, documentation |
|-------------------|--|

|                                 |                           |
|---------------------------------|---------------------------|
| Packaging Dimensions (LxWxH cm) | 33x45,6x10,2 (carry case) |
|---------------------------------|---------------------------|

|   |                  |
|---|------------------|
| TRIOS 6 scans are compatible with TRIOS Dx Plus | Yes <sup>2</sup> |
|---|------------------|

### HARDWARE

|        |           |
|--------|-----------|
| Weight | 308 grams |
|--------|-----------|

|                               |                  |
|-------------------------------|------------------|
| Scanner Dimensions (LxWxH cm) | 3.8 x 3.9 x 26.6 |
|-------------------------------|------------------|

|                           |                   |
|---------------------------|-------------------|
| Tip Dimensions (LxWxH cm) | 9,5 x 2,26 x 1,95 |
|---------------------------|-------------------|

### BATTERY

|                            |                   |
|----------------------------|-------------------|
| Battery Dimensions (LxWxH) | 7.7 x 2.17 x 2.42 |
|----------------------------|-------------------|

|                                      |   |
|--------------------------------------|---|
| Battery Duration (non-stop scanning) | Instant heating and up to 45 min of continuous scanning per battery |
|--------------------------------------|---|

|                       |  |
|-----------------------|--|
| Battery Charging Time | 60 minutes to achieve 80 % battery power. Less than 2 hours to achieve 100 % battery power |
|-----------------------|--|

<sup>1</sup> Tooth wear, gingival recession, plaque, surface and proximal caries. Proximal caries not AI detected.

<sup>2</sup> Please ask your reseller or 3Shape contact for availability of 3Shape products in your country. TRIOS Dx is not cleared by the FDA for clinical use in the USA.



Let's change dentistry together

3shape 