



What was once Hygenic is now HySolate the new Dental Dam Range from Coltene

Complete isolation with Dental Dam is part of the endodontic and restorative work flow. Teeth are isolated from the rest of the oral cavity and the working field is set dry. Using Dental Dam protects patients form aspiration or ingestion of foreign objects and protects soft tissue from contact with instruments.

HySolate Black Dental Dam has the template for marking the tooth position printed on. This saves the working step of marking the dam before punching. Like the Knight figure in a chess play you can leap steps, move faster and clever towards the goal of absolute moisture control.

The black colour provides ultimate, clear contrast and is ideal for taking photographs.

- Pre-printed dental dam reduces the number of working steps
- Black colour provides ultimate, clear contrast and is ideal for taking photographs.
- Powder free, low protein latex quality reduces the risk of developing latex hypersensitivity
- Acts as an infection control barrier
- Creates a dry operating field









FROM HYGENIC		TO HYSOLATE	
CW-H00539	HYGENIC Dental Dam Latex Dark Medium 152x152mm Box 36 sheets	CW-60031352	HySolate Dental Dam Black Latex Printed Medium 152x152 36 box
CW-H00540	HYGENIC Dental Dam Latex Dark Heavy152x152mm Box 36 sheets	CW-60031353	HySolate Dental Dam Black Latex Printed Heavy 152x152 36 box
CW-H00541	HYGENIC Dental Dam Latex Dark X Heavy152x152mm Box 36 sheets	CW-60031354	HySolate Dental Dam Black Latex Printed X-Heavy 152x152 36 box

1300 65 88 22 henryschein.com.au

HySolate Latex Dental Dam is made of pure, natural rubber latex and is powder free. Powder free, low protein Dental Dam is a simple and clever measure to reduce the risk of developing latex hypersensitivity. Latex Dental Dam provides strong retraction. The lighter colours have the advantage of naturally illuminating the operating field whereas the darker colours can help with visual contrast.

- Powder free, low protein latex dam reduces the risk of developing latex hypersensitivity
- Increased visibility on the treatment site through isolation
- Acts as an infection control barrier
- Creates a dry operating field



FROM HYGENIC		TO HYSOLATE	
CW-H00533	HYGENIC Dental Dam Latex Light Thin 152x152mm Box 36 sheets	CW-60034427	HySolate Dental Dam Light Ltx Thin 152x152 36 sheets
CW-H00534	HYGENIC Dental Dam Latex Light Medium 152x152mm Box 36 sheets	CW-60034428	HySolate Dental Dam Light Latex Medium 152x152 36 sheets
CW-H02146	HYGENIC Dental Dam Latex Green Thin 6" x6" Box of 36 sheets	CW-60034418	HySolate Dental Dam Green Latex Thin 152x152 36 sheets
CW-H02147	HYGENIC Dental Dam Latex Med Green 152x152mm Box 36 Sheets	CW-60034419	HySolate Dental Dam Green Latex Medium 152x152 36 sheets
CW-H11625	ELASTI-DAM Dental Dam PF Latex 152 x 152 Medium Box 36 Green	CW-60034419	HySolate Dental Dam Green Latex Medium 152x152 36 sheets
CW-H02148	HYGENIC Dental Dam Latex Green Heavy Box of 36 sheets	CW-60034420	HySolate Dental Dam Green Latex Heavy 152x152 36 sheets
CW-H11676	ELASTI-DAM Dental Dam PF Latex 152 x 152 Heavy Box 36 Green	CW-60034420	HySolate Dental Dam Green Latex Heavy 152x152 36 sheets
CW-H02149	HYGENIC Dental Dam Latex Hvy Green 152x152mm Box 36 Sheets	CW-60034421	HySolate Dental Dam Green Latex X-Heavy 152x152 36 sheets
CW-H03529	HYGENIC Dental Dam Latex Blue Colour Medium 127x127mm Box 52	CW-60034441	HySolate Dental Dam Blue Latex Medium 127x127 52 sheets
CW-H03530	HYGENIC Dental Dam Latex Blue Colour Medium 152x152mm Box 36	CW-60034442	HySolate Dental Dam Blue Latex Medium 152x152 36 sheets
CW-H11627	ELASTI-DAM Dental Dam PF Latex 152 x 152 Medium Box 36 Blue	CW-60034442	HySolate Dental Dam Blue Latex Medium 152x152 36 sheets
CW-H04247	HYGENIC Dental Dam Latex Green Colour Med 152x152mm Box 360	CW-60034423	HySolate Dental Dam Green Latex Medium 152x152 360 sheets
CW-H04639	FIESTA Dental Dam Latex Asstd Colours Thin 6"x 6" Box of 36	CW-60034433	HySolate Fiesta Dental Dam Latex Assorted Thin 152x152 36 sheets
CW-H04641	FIESTA Dental Dam Latex Asstd Colours Medium 5'x5" Box of 52	CW-60034435	HySolate Fiesta Dental Dam Latex Assorted Medium 127x127 52 sheets
CW-H04642	FIESTA Dental Dam Latex Asstd Colours Medium 6x6" Box of 36	CW-60034436	HySolate Fiesta Dental Dam Latex Assorted Medium 152x152 36 sheets
CW-H04644	FIESTA Dental Dam Latex Asstd Colours Heavy 6'x6" Box of 36	CW-60034438	HySolate Fiesta Dental Dam Latex Assorted Heavy 152x152 36 sheets
CW-H07315	HYGENIC Dental Dam Latex Blue Colour 152 x 152mm Box of 36	CW-60034443	HySolate Dental Dam Blue Latex Heavy 152x152 36 sheets
CW-H11624	ELASTI-DAM Dental Dam PF Latex 127x127mm Medium box 52 Green	CW-60034415	HySolate Dental Dam Green Latex Medium 127x127 52 sheets
CW-H11678	ELASTI-DAM Dental Dam PF Latex 152 x 152 Heavy Box 36 Blue	CW-60034415	HySolate Dental Dam Green Latex Medium 127x127 52 sheets



Scan here to view the HySolate Dental Dam Range on our website

1300 65 88 22 henryschein.com.au