




EdgeFile X7 GP & PP

Where Form Fitting Meets
Performance



Meet the new Gutta Percha Points and Absorbent Paper Points for Endodontic Treatment from EdgeEndo.

The new EdgeEndo X7 Gutta Percha range has been created to be form fitting, providing assurance when undertaking endodontic treatment after using the EdgeEndo X 7 range of NiTi rotary files.

The new EdgeEndo X7 Gutta Percha offers:

- Vivid length marks calibrated at 16mm, 19, 22mm and 24mm
- Requisite high level of radiopacity
- Accurate and form fitting sizes adhering to ISO standard sizing

All leading to a positive and predictable long term treatment outcome.

The new EdgeEndo absorbent paper points are conveniently packaged, and gamma ray sterilised for your peace of mind and infection control. They also offer vivid length marks calibrated at 16mm, 19, 22mm and 24mm and are lintless.



Code	Description
X7 EdgeFiles Taper 0.04	
EN-SX70421MM	Assorted
EN-SX7170421RF	Size 17
EN-SX7200421RF	Size 20
EN-SX7250421RF	Size 25
EN-SX7300421RF	Size 30
EN-SX7350421RF	Size 35
EN-SX7400421RF	Size 40
EN-SX7450421RF	Size 45
X7 EdgeFiles Taper 0.06	
EN-SX70621MM	Assorted
EN-SX7170621RF	Size 17
EN-SX7200621RF	Size 20
EN-SX7250621RF	Size 25
EN-SX7300621RF	Size 30
EN-SX7350621RF	Size 35
EN-SX7400621RF	Size 40
EN-SX7450621RF	Size 45

Code	Description
Gutta Percha 0.04 - 60 / pack	
EN-X7GP04AST-1	GP Assorted 0.04
EN- X71704GP-1	GP point 17.04
EN-X72004GP-1	GP point 20.04
EN-X72504GP-1	GP point 25.04
EN-X73004GP-1	GP point 30.04
EN-X73504GP-1	GP point 35.04
EN-X74004GP-1	GP point 40.04
EN-X74504GP-1	GP point 45.04
Gutta Percha 0.06 - 60 / pack	
EN-X7GP06AST-1	GP Assorted 0.06
EN-X71706GP-1	GP point 17.06
EN-X72006GP-1	GP point 20.06
EN-X72506GP-1	GP point 25.06
EN-X73006GP-1	GP point 30.06
EN-X73506GP-1	GP point 35.06
EN-X74006GP-1	GP point 40.06
EN-X74506GP-1	GP point 45.06

Code	Description
Paper Points 0.04 - 100 / pack	
EN-SPP04AST-1	PP Assorted 0.04
EN-SPP1504-1	PP 15.04
EN-SPP2004-1	PP 20.04
EN-SPP2504-1	PP 25.04
EN-SPP3004-1	PP 30.04
EN-SPP3504-1	PP 35.04
EN-SPP4004-1	PP 40.04
EN-SPP4504-1	PP 45.04
EN-SPP5004-1	PP 50.04
Paper Points 0.06 - 100 / pack	
EN-SPP06AST-1	PP Assorted 0.06
EN-SPP1506-1	PP 15.06
EN-SPP2006-1	PP 20.06
EN-SPP2506-1	PP 25.06
EN-SPP3006-1	PP 30.06
EN-SPP3506-1	PP 35.06
EN-SPP4006-1	PP 40.06
EN-SPP4506-1	PP 45.06
EN-SPP5006-1	PP 50.06