


# Essentia<sup>®</sup> Masking Liner

The Ultimate masking agent



## Covers amalgam staining

Essentia Masking Liner is an opaque, heavily-filled injectable composite that will block discolouration without increasing value, creating the base for great aesthetics.



Images courtesy of Prof M. Peumans

## Blocks discolouration

The Masking Liner will wet and adapt well to the base of cavities and is well suited to placement in very thin layers.



Images courtesy of Dr. Chee Chang

## Essentia<sup>®</sup> Masking Liner refill

1 x Syringe 3.4gm (2.0mL), 20x Dispensing Tips



Ordering code  
GC-900971

