

The Monthly Drill Down

Handpieces



BA Ultimate

German engineering.
Titanium strength.
Ultimate control.

Looking to step up from your everyday high-speed handpiece? The BA Ultimate range delivers powerful cutting performance with up to 22W of power, at a highly competitive price. With standard and mini-head options, it's a versatile solution designed to suit a wide range of clinical preferences and existing systems.

[Shop now](#)



Titanium body

Lightweight strength, made in Germany

Smart Coat finish

Scratch-resistant, better grip

Quadruple spray

Consistent cooling at the bur

2-year warranty

On BA Ultimate Power+

Our customers are loving



Morita TorqTech Blue Band >

Tech you can trust.

KaVo EXPERTmatic Blue Band >

Less effort. More smiles.

BA Blue Band >

Tough on work, easy on budget.

Practice trends we're seeing

Do you know how many handpieces your practice actually needs?

Reprocessing turnaround, rather than clinical wear, is increasingly determining how many handpieces a practice needs. Under AS 5369:2023, handpieces are unavailable for patient care during cleaning, lubrication and sterilisation. Practices are therefore sizing fleets around appointment volumes and reprocessing times, while planned servicing and loan units help minimise disruption.

Internal data indicates an average service turnaround of around seven days, and practices with at least one spare high-speed handpiece per chair, experience fewer scheduling disruptions. Servicing is as much a scheduling decision as a maintenance one.



Resources [Standards Australia, AS 5369:2023 - Reprocessing of reusable medical devices, 2023.](#)



Handpiece down? Hayes is up.

Factory-trained technicians, genuine OEM parts and small equipment rentals so your day doesn't stop.

[Book now >](#)

Simple maintenance tips for BA handpieces >

Product & servicing resources




The hidden danger in dental handpieces >

Maintenance mistakes that can lead to patient burns.

How to properly maintain and extend the life of KaVo instruments >

Benefits of automated handpiece lubrication with Morita's Lubrina 2 >

Clinical articles



Celebrating 20 years of Practice Pink worldwide



This year marks a major milestone for Practice Pink, Henry Schein's global initiative supporting cancer awareness, prevention, early detection, and patient care across the healthcare community.

More than \$3M AUD donated globally to help advance cancer awareness, prevention, early detection and patient care.

Over \$15,600 raised by our Practice Pink charity in 2025, for the Children's Cancer Foundation, thanks to the support of customers like you.



Support The Children's Cancer Foundation

On 7 - 8 October, we'll donate a percentage of profits from every eligible purchase + Spend \$1,000 get 5% off



Use code at checkout.

SUPPORTPINK

Excludes licensed product. Terms and conditions apply.

Where prevention meets peace of mind

Annual servicing isn't just a recommendation: it's a critical requirement to support operational safety, warranty policies, compliance, and long-term reliability.

Discover how Preventative Maintenance can add an extra layer of protection to your equipment investment — helping support long-term performance, reliability and peace of mind.

[View brochure >](#)



Future ready



3D Printing? We've got you covered.

More materials. More indications. More opportunities to grow.

From models and splints to denture bases, try-ins and permanent restorations, we offer one of Australia's most comprehensive 3D printing resin portfolios.

[Shop the range >](#)