

Pro Freehand Kit Drill Usage Chart

Description	Part Number	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Tapered Pilot Drill, 2.0mm	TSD2020PD																				
Tapered Dense Bone Drill, 2.5mm	TSD2025DB																				
Tapered Dense Bone Drill, 2.8mm	TSD2028DB																				
Tapered Soft Bone Drill, 2.8mm	TSD2028SB																				
Tapered Dense Bone Drill, 3.2mm	TSD2032DB																				
Tapered Soft Bone Drill, 3.2mm	TSD2032SB																				
Tapered Dense Bone Drill, 3.7mm	TSD2037DB																				
Tapered Soft Bone Drill, 3.7mm	TSD2037SB																				
Tapered Dense Bone Drill, 4.1mm	TSD2041DB																				
Tapered Soft Bone Drill, 4.1mm	TSD2041SB																				
Tapered Dense Bone Drill, 4.7mm	TSD2047DB																				
Tapered Soft Bone Drill, 4.7mm	TSD2047SB																				
Tapered Dense Bone Drill, 5.4mm	TSD2054DB																				
Tapered Short Pilot Drill, 2.2/3.2mm	TDS32PD																				
Tapered Short Soft Bone Drill, 2.4/3.31mm	TDS33SB																				
Tapered Short Dense Bone Drill, 3.7mm	TDS37DB																				
Tapered Short Crestal Bone Drill, 4.2mm	TDS42CB																				
Tapered Short Soft Bone Drill, 2.8/3.65mm	TDS37SB																				
Tapered Short Dense Bone Drill, 4.15mm	TDS42DB																				
Tapered Short Crestal Bone Drill, 4.6mm	TDS46CB																				
Tapered Short Soft Bone Drill, 3.2/4.16mm	TDS42SB																				
Tapered Short Dense Bone Drill, 4.8mm	TDS48DB																				
Tapered Short Crestal Bone Drill, 5.2mm	TDS52CB																				
Tapered Short Soft Bone Drill, 3.6/4.65mm	TDS47SB																				
Tapered Short Dense Bone Drill, 5.4mm	TDS54DB																				
Tapered Short Crestal Bone Drill, 5.8mm	TDS58CB																				

Important Note:
All drills should be replaced after they have been used 12 to 20 times, depending upon bone density.

Important Note: All drills should be replaced after they have been used 12 to 20 times, depending upon bone density.

customer care: 1300 769 861 • shop online at biohorizonscamlog.com.au

©BioHorizons Implant Systems Inc. All Rights Reserved.