


Tapered Pro Conical streamlined full-arch workflow

Fuel full-arch workflows with Tapered Pro Conical 3.8mm implant

The BioHorizons Tapered Pro Conical implant lineup offers the 3.8mm diameter implant in two prosthetic platforms: narrow and regular. This expanded portfolio streamlines full-arch workflows by allowing the use of a single prosthetic platform for any full-arch case.

surgical flexibility


Achieve optimized esthetic outcomes and conquer anatomical challenges by performing full-arch indications with all-yellow prosthetic connection implants or with all-gray prosthetic connection implants.



Narrow platform

simplified prosthetic workflow

All abutment components are interchangeable across all implant sites to reduce cost and inventory, minimize the chance of mismatched components, streamline lab communication, and enable faster turnaround during prosthesis fabrication.



Regular platform







streamlined surgical workflow

Reduce surgical chairtime by leveraging BioHorizons consolidated surgical kits designed for Tapered Pro Conical implants. Reduced instrumentation needs save space and cost.



Tapered Pro Conical

ordering information


						
Body diameter	3.3mm	3.8mm	3.8mm	4.2mm	4.6mm	5.2mm
Prosthetic connection	 Narrow	 Narrow	 Regular	 Regular	 Regular	 Regular
Platform switch	0.3mm	0.3mm	0.4mm	0.4mm	0.6mm	0.9mm
Laser-Lok® zone	1.8mm	1.8mm	1.8mm	1.8mm	1.8mm	1.8mm
Apical diameter	2.4mm	2.8mm	2.8mm	3.0mm	3.2mm	3.3mm
9.0mm length	CTA3309	CTA3809	CTA3809R	CTA4209	CTA4609	CTA5209
10.5mm length	CTA3310	CTA3810	CTA3810R	CTA4210	CTA4610	CTA5210
12.0mm length	CTA3312	CTA3812	CTA3812R	CTA4212	CTA4612	CTA5212
15.0mm length	CTA3315	CTA3815	CTA3815R	CTA4215	CTA4615	CTA5215
18.0mm length	CTA3318	CTA3818	CTA3818R	CTA4218	CTA4618	—



Pro Freehand Surgical Kit

FSK5000

Pro Freehand Surgical Kit



Guided Surgery Kit

PRO5000

Pro Guided Kit

BioHorizons® and Laser-Lok® are registered trademarks of BioHorizons Implant Systems Inc. Not all products are available in all countries. ©BioHorizons. All Rights Reserved.

BIOHORIZONS®



ML0348



REV A MAY 2025

customer care: 888.246.8338
shop online at store.biohorizons.com