

# Planmeca Compact™ iClassic



# The ideal choice for the needs of general dentistry

A cost-effective dental unit with all the essential functionalities

## Basic unit and instrument functions

*Planmeca Compact™ i Classic has a traditional keyboard with all the essential unit, instrument and infection control functions. A keyboard is also available for the assistant suction arm to provide easy and ergonomic use.*

### Durable aluminium console

The dental unit keyboard is integrated in the new and robust aluminium console. As surface material, aluminium is very durable and tolerates disinfectants used in the dental clinic well.



### Predefined instrument range


**Planmeca Compact™ i Classic** is offered with a predefined instrument range. Basic instrumentation covers all necessary instruments including integrated electric micromotors. All instrument adjustments can be done conveniently from the dental unit keyboard. The instrumentation can also be extended with electric scaler, LED polymerisation light and intraoral camera.

## Optimised ergonomics

*The functional design and carefully selected features of Planmeca Compact™ i Classic guarantee a smooth workflow and an excellent working position.*

### Perfectly balanced instrument arms

The balanced instrument arms guarantee ergonomic instrument use since only a tiny movement is sufficient to activate the instrument. Alternatively, Planmeca Compact i Classic can be equipped with hanging tube instruments.



### Functional design

The side-lifting chair and small cuspidor base maximise the legroom for both the dentist and assistant. The lean and narrow backrest design provides an excellent access to the treatment area.

### All operations from one foot control

All necessary functions can be carried out from the unit's robust foot control, keeping your hands free and minimising cross-contamination.



# Excellent patient comfort

The dental unit comes with a manual headrest and fixed legrest.  
 The double articulated headrest allows comfortable positioning of the patient's head.  
 Thanks to the smooth headrest length adjustment, the chair can be easily adapted for patients of all sizes. The chair with fixed legrest provides excellent patient comfort and uniform upholstery.



**Planmeca Compact™ i Classic** is available with both **Comfy™** and **Ultra Relax™** upholsteries in multiple colours.



*Planmeca Showroom*



*infokit.planmeca.com*



*Planmeca iRomexis™*



## PLANMECA

Asentajankatu 6 | 00880 Helsinki | Finland | tel. +358 20 7795 500 | fax +358 20 7795 555 | sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations.  
 Some products displayed above may not be available in all countries or areas. Rights for changes reserved.

Planmeca, All in one, Anatomat Plus, Cobra, Comfy, Digital perfection, Economat Plus, Elegant, Flexy, Perio Fresh, PlanEasyMill, Planmeca AINO, Planmeca ARA, Planmeca Chair, Planmeca Compact, Planmeca Intra, Planmeca iRomexis, Planmeca Lumion, Planmeca Minea, Planmeca Minendo, Planmeca Minetto, Planmeca Noma, Planmeca Online, Planmeca PlanCAD, Planmeca PlanMill, Planmeca PlanScan, Planmeca Planosil, Planmeca ProCeph, Planmeca ProFace, Planmeca ProlD, Planmeca ProMax, Planmeca ProModel, Planmeca ProOne, Planmeca ProScanner, Planmeca ProSensor, Planmeca ProX, Planmeca Romexis, Planmeca SingLED, Planmeca Sovereign, Planmeca Ultra Low Dose, Planmeca Vision, Planmeca Waterline Cleaning System, Proline Dental Stool, Saddle Stool, SmartPan, SmartTouch, Trendy, Ultra Relax and UltraTouch are registered or non-registered trademarks of Planmeca in various countries.

