

# Planmeca ProCeph

*One shot cephalostat*



*Planmeca ProCeph one shot cephalostat takes the image with one short exposure. It gives excellent diagnostic image quality of digital cephalometric radiographs with low dose.*

**PLANMECA** 



- For Planmeca ProMax, Planmeca ProMax 3D, Planmeca ProMax 3D s and Planmeca ProMax 3D Mid models
- Image sizes from 20 x 25 cm to 30 x 25 cm
- Excellent image quality with large, high-sensitivity flat-panel
- Exposure time 1 s – no motion artefacts
- Low dose
- Coming soon

Planmeca Oy  
Asentajankatu 6 | 00880 Helsinki | Finland  
tel. +358 20 7795 500 | fax +358 20 7795 555  
sales@planmeca.com | www.planmeca.com