



Durable polishers for direct and indirect restorations. Specifically designed to deliver superior high-gloss polish on composite, hybrid ceramics and zirconia.

Advantages

- Diamond-infused polishers that fulfill the highest industry standards for polishing materials
- Color coded instruments help identify the right instrument for the polish phases
- Preservation of the surface structure
- Less heat generation & excellent durability
- Designed for excellent polishing results in less time

Henry Schein
Brand



COPTER COMPOSITE POLISHERS - RA SHANK

- 2-step polishing system provides superior high-gloss finishes on composite and hybrid ceramic materials.
- Diamond impregnated, flexible spirals allow easy access to all tooth surfaces, including hard to reach interproximal areas.
- Red Pre-Polisher is designed to remove scratches and create initial shine, while the Blue High-polisher produces a final high-gloss polish.
- Available in 2 sizes - Large (14mm) for most general cases and Small (10mm) for more confined or hard-to-reach areas with limited access.
- Autoclavable.




10mm	10mm	14mm	14mm
Small	Small	Large	Large
Red Pre-Polish	Blue High-Polish	Red Pre-Polish	Blue High-Polish

Product	Quantity	Item Code	Contains
Starter Set	12 per Set	988-4956	3 of each: Small Pre-polish, Small High-polish, Large Pre-polish, Large High-Polish
Refill Set Small Size	6 per Pkg	988-4955	3 of each Small Pre-polish, Small High-Polish
Refill Set Large Size	6 per Pkg	988-4954	3 of each Large Pre-polish, Large High-Polish
Trial Set Small Size	2 per Pkg	988-4952	1 of each Small Pre-polish, Small High-Polish
Trial Set Large Size	2 per Pkg	988-4953	1 of each Large Pre-polish, Large High-Polish
Individual Size Refills	3 per Pkg	988-4957	Small Pre-Polish Red
Individual Size Refills	3 per Pkg	988-4958	Large Pre-Polish Red
Individual Size Refills	3 per Pkg	988-4959	Small High-Polish Blue
Individual Size Refills	3 per Pkg	988-4960	Large High-Polish Blue

COPTER COMPOSITE POLISHERS SINGLE USE -RA SHANK

- Single-use, 2-step polishing system for composite and hybrid ceramic materials.
- Diamond impregnated, flexible spirals allow easy access to all tooth surfaces, including hard to reach interproximal areas.
- Lime Green Pre-polisher is designed to remove scratches and create initial shine, while the Gray High-polisher produces a final high-gloss polish.
- Small (10mm) size allows easy access to all areas.
- Single-use allows for optimal performance with every patient.



10mm	10mm
Small	Small
Lime Green Pre-Polish	Gray High-Polish

Product	Quantity	Item Code	Contains
Starter Set	50 per Set	988-4969	25 of each: Pre-polish, High-Polish
Individual Size Refills	50 per Set	988-4970	Pre-Polish Lime Green
Individual Size Refills	50 per Set	988-4971	High-Polish Gray

COPTER ZIRCONIA POLISHERS - RA SHANK

- 2-step polishing system specifically designed for zirconia restorations.
- Highly concentrated diamond impregnated, provides ultimate efficiency and avoids repeated glazing after retouching.
- Orange Pre-polisher is designed to remove scratches and create initial shine, while the Dark Blue High-polisher produces a final high-gloss polish.
- Large (14mm) size is perfect for all cases.
- Used after grinding with Copter Zirconia Grinders.
- Autoclavable.



14mm	14mm
Large	Large
Orange	Dark Blue
Medium Pre-Polish	Fine Hih-Polish

Product	Quantity	Item Code	Contains
Starter Set	6 per Set	988-4963	3 of each: Pre-polish Medium, High-Polish Fine
Individual Size Refills	3 per Pkg	988-4964	Pre-Polish Orange
Individual Size Refills	3 per Pkg	988-4965	High-Polish Dark Blue

COPTER ZIRCONIA GRINDERS - RA SHANK

- Diamond impregnated, synthetic bonded instrument designed for adjusting zirconia restorations.
- Coarse grit for high abrasion with minimal heat development.
- 3 shapes provide versatility for all needs.
- Used prior to polishing with Copter Zirconia Polishers.
- Autoclavable.



Point	Cylinder	Disc
Grey	Grey	Grey

Product	Quantity	Item Code	Contains
Individual Sizes	3 per Pkg	988-4966	Point
Individual Sizes	3 per Pkg	988-4967	Cylinder
Individual Sizes	3 per Pkg	988-4968	Disc

COPTER ZIRCONIA ADJUSTMENT SET- RA SHANK

- Complete set for intra-oral adjustments and polishing of zirconia restorations.
- Autoclavable.

			14mm	14mm
Point	Cylinder	Disc	Large	Large
Grey	Grey	Grey	Orange	Dark Blue
			Medium Pre-Polish	Fine High-Polish

Product	Quantity	Item Code	Contains
Set	5 per Set	988-4962	1 of each: Copter Polishers- Pre-polish Medium, High-Polish Fine; Copter Grinders - Point, Cylinder, Disc.

COPTER ZIRCONIA ADJUSTMENT SET- HP SHANK

- Complete set for extra-oral adjustments and polishing of zirconia restorations.
- Autoclavable.

		14mm	14mm	14mm
Cylinder	Wheel	Large	Large	Large
Green	Green	Orange	Dark Blue	Dark Blue
Coarse	Coarse	Medium Pre-Polish	Fine High-Polish	Fine High-Polish

Product	Quantity	Item Code	Contains
Set	4 per Set	988-4961	1 of each: Copter Polishers - Pre-Polish Medium, High-Polish Fine; Copter Grinders- Wheel, Cylinder

Henry Schein
Brand

1300 65 88 22 | henryschein.com.au

SCAN CODE
FOR MORE DETAILS

