# **Sof-Lex**<sup>TM</sup>

#### Finishing and Polishing System

## Pick up the Pace without Sacrificing Great Results.

### **Indications**

- Sof-Lex polishing discs and finishing strips for polishing of composite, compomer, metal and ceramic.
- Sof-Lex polishing discs for facial and lingual margins of anterior and posterior restorations.
- Sof-Lex finishing strips for interproximal areas.

## **Advantages**

#### **Sof-Lex Polishing Discs:**

- Flexible polishing discs with aluminium oxide coating.
- Four different abrasive grades (coarse to superfine) and two disc thicknesses (normal and extra thin).
- Each abrasive grade is identified by a special colour coding system.
- Discs can be easily changed due to the patented pop-on mandrel system.
- Fracture resistant and sterilisable metal mandrel.



Sof-Lex Polishing Discs

#### **Sof-Lex Finishing Strips:**

- Thin flexible backing with aluminium oxide coating.
- Centre gapped strip with two abrasive grades per strip, which allows interproximal access.
- Three different types (fine/superfine, coarse/medium, coarse/medium (narrow)).

Abrasive Grades	Sof-Lex Moores Mandrel Type	Sof-Lex		Sof-Lex Extra Thir	n (XT)	Sof-Lex Finishing Strip
Superfine		•	0	•	•	SF
Fine	1958SF	1982SF	1981SF	2382SF	2381SF	F
	1958F	1982F	1981F	2382F	2381F	1956
Medium	1958M	1982M	1981M	2382M	2381M	M
Coarse	•	•	•	•	•	С
	1958C	1982C	1981C	2382C	2381C	1954N 1954